



BC-6800 Series Maintenance Report

RXQ-24032-BC-700 Series(V1.0)

Customer		Model	BC6800
Analyzer SN		Service Order	
Status	☐ Good		
Engineer		PM Date	
No.	Content	Acceptance	Result
1	Check the waste cistern floaters and WC1 drainage tube	WC1 & WC2 don't have much dirt on the floater.	☐ Done
2	Check the buffer bath	No crystals or liquid in the buffer bath.	☐ Done
3	Check tube pincher	Tube gripper is not deformed or bent.	☐ Done
4	Replace the blocking ring (if needed)	The blocking ring is not scratched or worn. Replace it if scratched or worn.	☐ Done
5	Check the drive pulley, driven pulley and slide rail	Drive pulley and Driven pulley can be rotated smoothly. Lubricate and clean if they are rusted.	☐ Done
6	Check the tube detection sensor status	Service-Debug & Self-test-Sensor Debug -Tube detection sensor -Inquire, the blank voltage rate is within 0.9~1.1.	☐ Done
7	Verify the sample probe position in the reaction bath	For the DIFF/RET bath, the probe shall fall in the center of the bath; but for the WNB bath, the probe should be about 2.5mm from the center of the bath.	☐ Done
8	Check the HGB and optical reaction bath bubble mixing tube	The HGB bath and WBC baths interior Teflon tube shall exceed the three-way connector for 4 mm. Their bubble mixing works normal and their waste can be discharged smoothly.	☐ Done
9	Check tube of reagent cap assembly	They are connected well without being bent or leakage.	☐ Done
10	Clean the crystal on the solenoid valve	No crystal or leakage on valves.	☐ Done
11	Clean the crystal on the syringe	No crystal or leakage on syringes.	☐ Done
12	Maintain the pinch valve tube	No leakage or damage on the pinch valve tubes.	☐ Done
13	Clean the piercing probe, probe wipe and probe wipe barrier	Probe & Probe wipe & Probe wipe barrier are clean and dry.	☐ Done

14	Replace the diluent filter	Replace diluent filter.	☐ Done
15	Check piercing probe abrasion status	The piercing probe is not worn out.	☐ Done
16	Clean the HGB bath RBC bath and optical reaction bath	HGB bath /RBC bath / WBC baths covers are clean.	☐ Done
17	Lubricate the moving parts	Mechanical moving parts can be move smoothly.	☐ Done
18	Check and replace the wearing parts	☐ Piercing probe and probe wipe ☐ Pinch valve tubes from PV5 & PV6 ☐ Probe wipe filter ☐ Waste filter ☐ Sheath fluid filter ☐ Air pump	☐ Done
19	Perform Probe Cleanser maintenance to the waste cistern (if necessary)	WC1 & WC2 don't have much dirt on the floater.	☐ Done
Overall conclusion of installation acceptance		☐ Pass	