



Installation Checklist

Work Order		Ppm Date	
Analyzer Model		Hospital	
Serial Number		Department	
Windows OS		SW Version	

No. NO.	Items Inspection Item		Description Description	Reference Reference	Result Result	Re- marks
1	Operating environment	Ambient temperature	Within the normal operating temperature range	10°C□35°C	□Pass □Fail	
2		Humidity	Within the normal operating humidity range	30%□85%	□Pass □Fail	
3		Atmospheric pressure	Within the normal operating atmospheric pressure range	70KPa□106KPa	□Pass □Fail	
4		Electromagnetic interference	Away from electromagnetic interference source	/	□Pass □Fail	
5	Space Requirements		Space	Meet the space requirements	□Pass □Fail	

No. NO.	Items Inspection Item	Description Description	Reference Reference	Result Result	Re- marks
6	Heat dissipation	Space requirement for radiator fan	Distance to wall: $\geq 250\text{mm}$	☐ Pass ☐ Fail	
7	Power requirements of the main unit	Power supply voltage	100V-240V~ (allowance $\pm 10\%$)	☐ Pass ☐ Fail	
8		Zero-to-earth voltage	< 5v	☐ Pass ☐ Fail	

No. NO.	Items Inspection Item	Description Description	Reference Reference	Result Result	Re- marks
1.	The package box and appearance of the main unit	Package and appearance	Package and appearance are in good conditions	☐ Pass ☐ Fail	
2.	Package and appearance of the autoloader	Package and appearance	Package and appearance are in good conditions	☐ Pass ☐ Fail	
3.	Packing list	Items in packing list	All items included	☐ Pass ☐ Fail	
4.	Removal of fixing parts for analyzer transport	Sampling assembly and blending assembly	Cut the wire ties	☐ Pass ☐ Fail	
5.		Fluorescent reagent container	Cut the wire ties	☐ Pass ☐ Fail	
6.		Rotary Compressing Assembly	Cut the wire ties	☐ Pass ☐ Fail	
7.		Sample Compartment	Cut the wire ties	☐ Pass ☐ Fail	
8.		Peripheral Blood Blending Assembly	Cut the wire ties	☐ Pass ☐ Fail	
9.	HGB blank voltage	HGB background voltage values	4.3-4.5 V	☐ Pass ☐ Fail	

No. NO.	Items Inspection Item	Description Description	Reference Reference	Value Vallue	Result Result	Re- marks
1.	Calibration factors	WBC	75□125		☐ Pass ☐ Fail	
		RBC	75□125		☐ Pass ☐ Fail	
		HGB	75□125		☐ Pass ☐ Fail	
		MCV	95□105		☐ Pass ☐ Fail	
		PLT	75□125		☐ Pass ☐ Fail	
		RBC-O Whether to run ☐ Yes ☐ No	75□125		☐ Pass ☐ Fail	
		PLT-O Whether to run ☐ Yes ☐ No	75□125		☐ Pass ☐ Fail	
		PLT-H Whether to run ☐ Yes ☐ No	75□125		☐ Pass ☐ Fail	
2.	Background (OV-WB)	WBC	$\leq 0.1 \times 10^9 / L$		☐ Pass ☐ Fail	
		RBC	$\leq 0.02 \times 10^{12} / L$		☐ Pass ☐ Fail	
		HGB	$\leq 1 g / L$		☐ Pass ☐ Fail	
		PLT	$\leq 3 \times 10^9 / L$		☐ Pass ☐ Fail	
3.	Optical gain confirmation	DIFF-FS	$\leq 5.0\%$		☐ Pass ☐ Fail	
		DIFF-SS	$\leq 5.0\%$		☐ Pass ☐ Fail	
		DIFF-FL	$\leq 5.0\%$		☐ Pass ☐ Fail	
		RET-FS	$\leq 10.0\%$		☐ Pass ☐ Fail	

No. NO.	Items Inspection Item	Description Description	Reference Reference	Value Vallue	Result Result	Re- marks
		Whether to run ☐ Yes ☐ No				
		RET-SS Whether to run ☐ Yes ☐ No	≤10.0%		☐ Pass ☐ Fail	
		RET-FL Whether to run ☐ Yes ☐ No	≤10.0%		☐ Pass ☐ Fail	
4.	Carryover in blood routine	WBC	≤1.0□		☐ Pass ☐ Fail	
		RBC	≤1.0□		☐ Pass ☐ Fail	
		HGB	≤1.0□		☐ Pass ☐ Fail	
		HCT	≤1.0□		☐ Pass ☐ Fail	
		PLT	≤1.0□		☐ Pass ☐ Fail	
5.	Repeatability (OV-WB)	WBC	≤2.5%		☐ Pass ☐ Fail	
		RBC	≤1.5%		☐ Pass ☐ Fail	
		HGB	≤1.0%		☐ Pass ☐ Fail	
		MCV	≤1.0%		☐ Pass ☐ Fail	
		HCT	≤1.5%		☐ Pass ☐ Fail	
		PLT	≤4.0%		☐ Pass ☐ Fail	
		ESR	≤1.0(SD)		☐ Pass ☐ Fail	
		RET# Whether to run	≤15%		☐ Pass ☐ Fail	

No. NO.	Items Inspection Item	Description Description	Reference Reference	Value Vallue	Result Result	Re- marks
		☐ Yes ☐ No				
		RET% Whether to run ☐ Yes ☐ No	≤15%		☐ Pass ☐ Fail	
6.	Repeatability (OV-PD) Whether to run ☐ Yes ☐ No	WBC	≤3.5%		☐ Pass ☐ Fail	
		RBC	≤2.0%		☐ Pass ☐ Fail	
		HGB	≤2.0%		☐ Pass ☐ Fail	
		MCV	≤3.0%		☐ Pass ☐ Fail	
		HCT	≤3.0%		☐ Pass ☐ Fail	
		PLT	≤8.0%		☐ Pass ☐ Fail	
		RET# Whether to run ☐ Yes ☐ No	≤30.0%		☐ Pass ☐ Fail	
		RET% Whether to run ☐ Yes ☐ No	≤30.0%		☐ Pass ☐ Fail	
7.	Background (OV-BF) Whether to run ☐ Yes ☐ No	WBC-BF	≤0.001×10 ⁹ /L		☐ Pass ☐ Fail	
		RBC-BF	≤0.003×10 ¹² /L		☐ Pass ☐ Fail	
8.	Carryover (OV-BF) Whether to run ☐ Yes ☐ No	WBC-BF	≤0.3%		☐ Pass ☐ Fail	
		RBC-BF	≤0.3%		☐ Pass ☐ Fail	
9.	Repeatability (OV-BF) Whether to run ☐ Yes ☐ No	WBC-BF	≤30%		☐ Pass ☐ Fail	
		RBC-BF	≤40% or ≤7000/μL		☐ Pass ☐ Fail	