

ERYTRA EFLEXIS

PREVENTIVE MAINTENANCE

Form version 1.0

40.000 tests
6 Years

Instrument:	Form Version: 1.0
Serial Number: EF-12345	Code: EFX-6Y
Installation Date: 2024-03-10	Maintenance Date: 2026-10-04

Customer Site: King Fahad Hosp.	Customer Name: Hematology Lab
Department:	
Tel:	E-Mail:

Service Engineer Name: A. Alstath	
Department:	
Tel:	E-Mail:

Number of tests of last year PM: 38120
Current number of tests: 40210

PARTS REPLACEMENT			Instrument OFF
ACTION	Checked/Done	Value	Comment
Probe (PR)	☒		
Toothed belt T2.5/245 (PIP Y-axis)	☒		
Toothed belt T2.5/245 (PIP X-axis)	☒		
Toothed belt T2.5/245 (CAT Y-axis)	☒		
Toothed belt T2.5/245 (CAT X-axis)	☒		
Aspiration pump VPB	☒		
VISUAL / PHYSICAL INSPECTIONS & CLEANING			Instrument OFF
ACTION	Checked/Done	Value	Comment
GENERAL			
General appearance	☒		
Output cables and connectors	☒		
Grounding continuity test. Switch off and unplug the electrical appliance. Check with a multimeter electric continuity between earth terminal of IEC inlet appliance and:			
Screw attaching centrifuge 1 to chassis	☒		
Screw attaching centrifuge 2 to chassis	☒		
MONITOR ARM			
Check that the arm movement is smooth	☒		
Check that the arm blocking screw prevents undesired impacts of the monitor against the instrument chassis	☒		
DOOR			
Check that movement is smooth	☒		
Check that door does not hit during automatic opening	☒		
Check that the latch enters smoothly when closing	☒		
OTHER			
If UPS is connected, check battery level	☒		
CARDS MODULE GATE			
Check that movement is smooth	☒		
SAMPLES AND REAGENTS ENTRY (SRE)			
Check sample tube holders	☒		
Check reagent vial/bottle lodgings	☒		
CARD DRAWERS			
Check card racks centerers and retainers	☒		
DILUTION CUP			
Check dilution cup drain tubing	☒		
CENTRIFUGE MODULE			
Removing centrifuge pot			

Clean* pot of the centrifuge	✓		
Check rotor and card holders status	✓		
Check that rotor movement is smooth	✓		
PIP MODULE			
Check belts tension and status	✓		
Check chains, wiring and tubing status	✓		
Check/clean* barcode reader mirrors	✓		
CAT MODULE			
Check belts tension and status	✓		
Check chains and wiring status	✓		
Check gripper pins, centerer and fingers movement/status	✓		
Check/clean* barcode reader mirror	✓		
Check/clean* presence detector	✓		
FLUIDICS MODULE			
Check status of the fluidic input/output system for each bottle pos.	✓		
Check status of the bottles	✓		
Perform a fluidic circuit general cleaning*	✓		
Check tubing	✓		
Check the solution lines filters status	✓		
Check vacuum filters status	✓		
Check that the main fluidic manifold has no scratches nor leakages	✓		
READER MODULE			
Removing reader cover			
Check/clean* reader pattern	✓		
Check/clean* reader camera optics	✓		
Check reader foam status	✓		
*If cleaning is required, please note it in the comments box.			
LUBRICATIONS			Instrument OFF
ACTION	Checked/Done	Value	Comment
DOOR			
Lubricating general module			
Lubricate door spring and door guides (grease)	✓		
CARD DRAWERS			
Lubricating general module			
Lubricate drawers latch (grease)	✓		
Lubricate linear guides (grease)	✓		
PIP MODULE			
Lubricating pipetting arm module			
Lubricate X-axis linear guide (oil)	✓		
Lubricate Y-axis linear guide (oil)	✓		
Lubricate homs linear guide (oil)	✓		
Lubricate probe leadscrew (grease)	✓		
Lubricate probe washer holder cam (grease)	✓		
Lubricate homs rack and pinion system (grease)	✓		
CAT MODULE			
Lubricating card transport arm module			
Lubricate X-axis linear guide (oil)	✓		
Lubricate Y-axis linear guide (oil)	✓		
Lubricate gripper linear guide (oil)	✓		
Lubricate gripper leadscrew (grease)	✓		
Lubricate gripper pins (oil)	✓		
Lubricate gripper centering system	✓		
FLUIDICS MODULE			
Lubricating fluidics module			
Lubricate bottle pin entry for each bottle position (grease)	✓		
GENERAL			
Check the general appearance of the instrument and if necessary clean the affected parts.	✓		
FUNCTIONAL CHECKINGS & ADJUSTMENTS			Instrument ON
ACTION	Checked/Done	Value	Comment
DILUTION CUP			

[Check proper stirrer operation](#)
(PROFILE 0: 400 ± 40 rpm)



CENTRIFUGE MODULE

[Centrifugation verification \(PROFILE 0:](#)
[1206 rpm ± 36 rpm with duration](#)
[of 550 s ± 10 s\)](#)



PIP MODULE

[Check presence detector alignment](#)



[Verify the CPD sensor with Leak-](#)
[tightness test in TSS](#)



[Check probe offset in TSS](#)



CAT MODULE

[Check gripper-presence detector](#)
[alignment](#)



READER MODULE

[Removing reader cover](#)

[Perform reader adjustment](#)



[Perform gripper adjustment in card](#)
[reader area](#)



SOFTWARE

[Folders inspection](#)



COMPUTER AND TOUCHSCREEN

[SSD available space](#)



[Touchscreen operability \(touch](#)
[functionality in all screen areas\)](#)



QUALIFICATION

ACTION	Checked/Done	Value	Comment
Perform Operational Qualification (DG- AT-TSM-000150)			