

ERYTRA EFLEXIS

PREVENTIVE MAINTENANCE

Form version 1.0

40.000 tests
6 Years

Instrument:	Form Version: 1.0
Serial Number:	Code:
Installation Date:	Maintenance Date:

Customer Site:	Customer Name:
Department:	
Tel:	E-Mail:

Service Engineer Name:	
Department:	
Tel:	E-Mail:

Number of tests of last year PM:
Current number of tests:

PARTS REPLACEMENT Instrument OFF

ACTION	Checked/Done	Value	Comment
Probe (PR)	☐		
Toothed belt T2.5/245 (PIP Y-axis)	☐		
Toothed belt T2.5/245 (PIP X-axis)	☐		
Toothed belt T2.5/245 (CAT Y-axis)	☐		
Toothed belt T2.5/245 (CAT X-axis)	☐		
Aspiration pump VPB	☐		

VISUAL / PHYSICAL INSPECTIONS & CLEANING Instrument OFF

ACTION	Checked/Done	Value	Comment
GENERAL			
General appearance	☐		
Output cables and connectors	☐		
Grounding continuity test. Switch off and unplug the electrical appliance. Check with a multimeter electric continuity between earth terminal of IEC inlet appliance and:			
Screw attaching centrifuge 1 to chassis	☐		
Screw attaching centrifuge 2 to chassis	☐		
MONITOR ARM			
Check that the arm movement is smooth	☐		
Check that the arm blocking screw prevents undesired impacts of the monitor against the instrument chassis	☐		
DOOR			
Check that movement is smooth	☐		
Check that door does not hit during automatic opening	☐		
Check that the latch enters smoothly when closing	☐		
OTHER			
If UPS is connected, check battery level	☐		
CARDS MODULE GATE			
Check that movement is smooth	☐		
SAMPLES AND REAGENTS ENTRY (SRE)			
Check sample tube holders	☐		
Check reagent vial/bottle lodgings	☐		
CARD DRAWERS			
Check card racks centerers and retainers	☐		
DILUTION CUP			
Check dilution cup drain tubing	☐		
CENTRIFUGE MODULE			
Removing centrifuge pot			

Clean* pot of the centrifuge	☐
Check rotor and card holders status	☐
Check that rotor movement is smooth	☐
PIP MODULE	
Check belts tension and status	☐
Check chains, wiring and tubing status	☐
Check/clean* barcode reader mirrors	☐
CAT MODULE	
Check belts tension and status	☐
Check chains and wiring status	☐
Check gripper pins, centerer and fingers movement/status	☐
Check/clean* barcode reader mirror	☐
Check/clean* presence detector	☐
FLUIDICS MODULE	
Check status of the fluidic input/output system for each bottle pos.	☐
Check status of the bottles	☐
Perform a fluidic circuit general cleaning*	☐
Check tubing	☐
Check the solution lines filters status	☐
Check vacuum filters status	☐
Check that the main fluidic manifold has no scratches nor leakages	☐
READER MODULE	
Removing reader cover	
Check/clean* reader pattern	☐
Check/clean* reader camera optics	☐
Check reader foam status	☐
*If cleaning is required, please note it in the comments box.	
LUBRICATIONS	
ACTION	Checked/Done
DOOR	
Lubricating general module	
Lubricate door spring and door guides (grease)	☐
CARD DRAWERS	
Lubricating general module	
Lubricate drawers latch (grease)	☐
Lubricate linear guides (grease)	☐
PIP MODULE	
Lubricating pipetting arm module	
Lubricate X-axis linear guide (oil)	☐
Lubricate Y-axis linear guide (oil)	☐
Lubricate homs linear guide (oil)	☐
Lubricate probe leadscrew (grease)	☐
Lubricate probe washer holder cam (grease)	☐
Lubricate homs rack and pinion system (grease)	☐
CAT MODULE	
Lubricating card transport arm module	
Lubricate X-axis linear guide (oil)	☐
Lubricate Y-axis linear guide (oil)	☐
Lubricate gripper linear guide (oil)	☐
Lubricate gripper leadscrew (grease)	☐
Lubricate gripper pins (oil)	☐
Lubricate gripper centering system	☐
FLUIDICS MODULE	
Lubricating fluidics module	
Lubricate bottle pin entry for each bottle position (grease)	☐
GENERAL	
Check the general appearance of the instrument and if necessary clean the affected parts.	☐
FUNCTIONAL CHECKINGS & ADJUSTMENTS	
ACTION	Checked/Done
DILUTION CUP	

[Check proper stirrer operation](#)
(PROFILE 0: 400 ± 40 rpm)



CENTRIFUGE MODULE

[Centrifugation verification \(PROFILE 0:](#)
[1206 rpm ± 36 rpm with duration](#)
[of 550 s ± 10 s\)](#)



PIP MODULE

[Check presence detector alignment](#)



[Verify the CPD sensor with Leak-](#)
[tightness test in TSS](#)



[Check probe offset in TSS](#)



CAT MODULE

[Check gripper-presence detector](#)
[alignment](#)



READER MODULE

[Removing reader cover](#)

[Perform reader adjustment](#)



[Perform gripper adjustment in card](#)
[reader area](#)



SOFTWARE

[Folders inspection](#)



COMPUTER AND TOUCHSCREEN

[SSD available space](#)



[Touchscreen operability \(touch](#)
[functionality in all screen areas\)](#)



QUALIFICATION

ACTION	Checked/Done	Value	Comment
Perform Operational Qualification (DG- AT-TSM-000150)	☐		