

ACL AcuStar/BIO-FLASH Installation Checklist

Account:		Date:	
Address:			
Serial No:			



NOTE: Do not perform any alignments until **step 20**.

Step	Action to Perform	Done
1.	Check the accessories for completeness.	
2.	Check the instrument, computer, touch screen, keyboard, and mouse for damage.	
3.	Remove the instrument and electronics enclosure covers.	
4.	Check that the instrument Controller and Motor Driver PCBs are properly seated.	
5.	Reinstall the electronic enclosure cover.	
6.	Connect the computer to the instrument (serial cable).	
7.	Install the Condensation Drainage kit (included in the box). If the kit is NOT GOING TO BE INSTALLED, then you must push the "Y" connector below the Reagent Carousel to the rear lower compartment. See Technical Note TN-HW-0131107B-Condensation drainage kit included in the spare box.	
8.	Power ON the Instrument and Computer. Do not initialize the instrument at this point!	
9.	When the Operation Software starts, log in as Service (fservice, or ilservice).	
10.	Check the touch screen, keyboard, and mouse for proper operation.	
11.	Check the AcuStar/BIO-FLASH software version (→ About).	
12.	Check the Instrument serial number (→ About).	
13.	Load the consumables (rinse, triggers, cuvettes).	
14.	Make sure the Cuvette waste bin is in the Cuvette waste drawer.	
15.	Initialize and Prime the instrument (x4).	
16.	Clean the Luminometer injector block. See Cleaning the Trigger Heater Assembly Injector Block on page 694.	
17.	Check all fluid connections for leaks.	
18.	Close the AcuStar/BIO-FLASH software.	
19.	Open the Field Service software (FSSW) (ECDiagnostics.exe file.)	

Step	Action to Perform	Done
20.	1. Reagent Pipettor to Incubator/Distributor. 2. Sample Pipettor to Incubator Distributor 3. Reagent Pipettor to Reagent Carousel 4. Sample Pipettor to Reagent Carousel 5. Sample Pipettor to Sample Carousel 6. MagWash Aspirator 7. MagWash LLS Calibration Perform the following alignments. (Keep the order shown.)	
21.	Perform a Flash backup. See Saving Flash Memory Contents on page 340.	
22.	Exit FSSW and run the AcuStar/BIO-FLASH Software.	
23.	Log in as "Manager" or "Specialist."	
24.	Install the required assay libraries (Assay Definitions→ Import).	
25.	Check the temperatures using the Status screen. • Incubator (36.0°C to 38.0°C) • Reagent (2.0°C to 6.0°C) • Trigger (39.0°C to 41.0°C)	
26.	Re-install the instrument covers.	
27.	Initialize and Prime the instrument.	
28.	Use customer labels to check the Sample Barcode Reader.	
29.	Run Verification Reagent Calibration	
30.	Run Verification Reagent QC	
31.	Qualify the instrument according to the Verification spreadsheet. ¹	

Customer Name:		Signature:	
Service Engineer Name:		Signature:	

¹Contact the Instrument provider for the current Verifi