



LIAISON XL (DiaSorin)

Fully Automated Chemiluminescence Analyser

CHECK LIST FOR DE-INSTALLATION



Abdulla Fouad for Medical

Supplies and Service

WO No: _____ Service Date: _____

Customer Name: _____ Department: _____

Instrument Serial No./ AFCO ID: _____ Instrument SW version: _____

UN-INSTALLATION PROCEDURE

No.	Procedure	Flag	Remark
1	Data deleting (see Chapter 4.3)		
1.1	Delete patient database.		
2	Decontamination (see Chapter 4.3)		
2.1	Empty Liquid Waste Containers, add hypochlorite solution (refer to chapter 5 of LIAISON XL User Manual) and load them back into the system.		
2.2	Run Service Decontamination2 task from the dedicated window of Main User Interface (Main category Maintenance, sub-category Schedule).		
2.3	Switch the liquid waste tank.		
2.4	Run Clean-up Washer task from the dedicated window of Main User Interface (Main category Maintenance, sub-category Schedule), in order to decontaminate the second liquid waste line also.		
2.5	Empty starter, cleaning, water and wash containers and intermediate containers.		
2.6	Load empty containers into the system and empty liquid lines (using the procedures of chapters 3.3.8.1, 3.3.8.2, 3.3.8.3 and 3.3.8.4).		
2.7	Remove all containers from the system.		
2.8	If not disposed, decontaminate liquid waste containers (see Chapter 4.3.1).		
2.9	Remove all cuvettes from the system (see Chapter 4.3.2).		
2.10	Remove all tubes out of the sample racks.		
2.11	Remove all integrals out of the reagent module.		
2.12	Remove all ancillaries out of the ancillary bay.		
2.13	Remove all disposable tips out of the disposable tip drawers.		
2.14	Remove and dispose the solid waste bag.		
2.15	Remove and dispose (if possible) the waste tank bulk (if previously installed).		
2.16	Remove and dispose (if possible) the solid waste bin (if previously installed).		
2.17	Remove and dispose the cuvette waste tube.		

2.18	Remove and dispose the wash/ waste station.		
2.19	Clean the pipettor, the analyzer interior, the primary and intermediate tanks and the covers (refer to chapter 8 of LIAISON XL User Manual)		
2.20	Perform additional cleanings of contaminated parts (see Chapter 4.3.4).		
2.21	Perform additional cleanings of not contaminated parts (see Chapter 4.3.4).		
2.22	Complete the certificate of the decontamination (included in DiaSorin Common Technical Note 1).		
2.23	Add the certificate of decontamination to the transportation documents.		
2.24	Make sure Reagent Area & Incubator Module top covers are available with instrument		

Page 1 | 2

PACKING

No.	Procedure	Flag	Remark
3	System Packing (see Chapter 4.4 for LIAISON XL and Chapter 25.5.3/ Chapter 25.5.4 for LIAISON XL LAS)		
3.1	Prepare the system for the packing.		
3.2	Pack the housing parts.		
3.3	Pack the covers.		
3.4	Pack the accessories as per the packing list		
3.5	Pack the UPS		
3.6	Pack the Printer		
3.7	Pack the LIS connectivity module/DiGi box if applicable		
4	Comment / Observation:		

Engineer Name: _____ Signature: _____ Date: _____