

CL-1200i Installation and Validation Report

1. Instrumental Qualification:

1.1 Dealer and Customer Information:

Dealer		Customer	
Name:	AFMS	Name:	
Contact Person:		Contact person:	
Telephone No:		Telephone No:	

1.2 Instrument Information:

Type of Instrument:			
Serial No.:		Installation Date:	

1.3 Transport Damage:

Package Damaged:	☐ Yes	☒ No
Machine (Or Accessory) Damaged:	☐ Yes	☒ No
If yes, please specify:		
Please return the damaged parts, including a malfunction report.		

1.4 Initial Check:

Parts missing? (Check against the packing list)			
☒ No		☐ Yes→Please specify:	
Part number:		Part description:	
Part number:		Part description:	
Part number:		Part description:	
Part number:		Part description:	

1.5 Installation:

Perform the next steps:
☒ Install the PC (<i>Computer, Monitor, Printer, Keyboard and Mouse</i>)
☒ Install the Software (<i>the software version: (V 01.02.33)</i>)
☒ Connect the PC, Analyzer, Monitor, <i>Keyboard and Mouse</i> .
☒ Remove the fixed plate.
☒ Install the Probe
☒ Connect the Fluidic Tubing and Prepare the Waste Liquid Tank
☒ Load the Wash Buffer and the Waste Cuvette Boxes
☒ Prepare the Cuvette Boxes and the Substrate Tubing Wash Solution

2. Operational Qualification:

Serial Number:			
☒ Power Up PC, and Exit to Desktop.			
☒ Power Up the analyzer, and run software from start menu.			
☒ Set the Atmosphere Pressure and Run "Fluidic Prime".			
☒ Confirm the Mechanical Position			
☒ Set the Initial Maintenance Record			
☒ Load the Reagents and Import the Items Parameters			
☒ Confirm the Basic Performance Testing			
Remarks:			
Technical acceptance: accepted			
Date:			
Service Engineer		Customer	
Name:		Name:	

3. Performance Qualification:

Clinical Performance Validation				
☒ Run Calibrations				
☒ Run QC				
Test Items	QC range level 1	Test result	QC range level 2	Test result
Remarks:				
Technical acceptance				
Date:				
Service Engineer			Customer	
Name:		Name:		