



ETI MAX 3000 (DiaSorin)
Fully Automated ELISA Analyser
Check list for New Installation



Report No: _____ Service Date: _____

Customer Name: _____ Department: _____

Instrument Serial No: _____ Instrument SW version: _____

INSTALLATION QUALIFICATION

No.	Procedure	Flag	Remark
1	Unpacking		
1.1	Unpack the complete system		
2	Inventory		
2.1	Check accessories according to the packing list		The content of the boxes should match with the description available on the Service Manual, chapter 3.4 on page 3-2
3	Installation		
3.1	Remove left and right rubber foam protections		
3.2	Remove plate transport protection		
3.3	Remove pipettor protections		
3.4	Remove rubber foam from wash buffer bottles area		
3.5	Remove washer head protection		
3.6	Remove reader plate drive protection		
3.7	Install the tip release panel and waste bag		
3.8	Install washer head manifold		
3.9	Install the Z axis shaft and the guiding curve		
3.10	Install the tip adapter and set the correct pipettor tubing length		
3.11	Install and connect the system liquid container		
3.12	Install and connect the waste container		
3.13	Install side covers and top cover		
3.14	Plug in instrument, PC and printer		
3.15	Connect the PC to the instrument (COM 1)		
3.16	Install the UPS if present		
3.17	Check that all the electrical connections are compliant with the applicable safety specifications		
3.18	Perform a general inspection on the main mechanical parts		
4	Preliminary Check		
4.1	Switch on instrument, PC and printer		
4.2	Perform a self-test for all the modules		

OPERATIONAL QUALIFICATION

No.	Procedure	Flag	Remark
5	Module Test		
5.1	Run the COP module test		
5.2	Run the rack area module test		
5.3	Run the reader module test		Including alignment plate, black plate and verification plate run
5.4	Run the washer module test		Including alignment plate and calibration jig run
5.5	Check that the vacuum sensor is correctly set		Use internal waste bottles to test it
5.6	Run the incubator module test		Including the shake mechanism test
5.7	Through the service software find out the shortest plate carrier		
5.8	Check plate transport teaching to detect any bad adjustment		Use the put and get or cycle features to test the positions. Make corrections to the teaching if required.
5.9	Run the plate transport module test		Including P test
6	Pipettor Module Test		
6.1	Preparation for pipettor teaching: - Inspect all racks (spring and barcode labels) - Check the correct 4PSCoo.INI setting - Start the teacher software		Restart the Teacher software if changes are applied to the 4PSCoo.INI file
6.2	Check the main reference position (Ref. Left)		
6.3	Check the tip eject position		
6.4	Check the wash/waste station positions		
6.5	Check the five tip racks reference positions		
6.6	Check the tip type determination position		
6.7	Check the rack area reference positions		
6.8	Check the four pipette areas reference positions		
6.9	Check the dilution area reference position		Use dilution tubes (12x75) with 100 μ l as dead volume for Zmax
6.10	Check the 6 bottles positions		
6.11	Check the correct pipettor alignment on sample and dilution tubes through the Service Software > Pipettor > LLD Test		Check then Rack area by three sample racks loaded in lanes 1, 6 and 12. Check the dilution area by dilution tubes loaded in positions 1 to 4, 9 to 12, 17 to 20, 141 to 144, 149 to 152 and 157 to 160.
6.12	Run the pipettor automatic module test		Including "Prime", "Endurance run" and "Pick-up and Eject Tips". Do not run the LLD-Test.
7	Washer Adjustment		
7.1	Check the washer default heights adjustment from the System Set-Up > Washer menu in the user software.		
7.2	Check that the default heights adjustments are suitable for all plates in use on the system		Particular adjustments for different types of plates can be applied in the assay file, to the washing step.

8	Final Tests		
8.1	QC kit run preparation: - Rename the Reagent.db file in use as Reagent.org - Copy the Reagent.db file from the mini CD supplied with the QC kit into the System subfolder - Copy the four assay files from the mini CD supplied with the QC kit into the assay subfolder		Content of the QC kit mini CD: Reagent.db, 8x100.asy, 100&50.asy, Dilution.asy, Val-Wash.asy
8.2	Run the four assays (8x100; 100&50; Dilution; Val-Wash) according to the instructions reported in the kit IFU		
8.3	Remove all files used only for PM verification purposes		
8.4	Restore the original Reagent.dbfile		
8.5	Back-up all calibrations and global files		*.dat; *.rac; *.mpc; *.glb

DOCUMENTATION/ REPORTS

No.	Procedure	Flag	Remark
9	Attachments		
9.1	Attach all module test reports: - COP - Rack area - Reader - Washer - Incubator - Plate transport - Pipettor		
9.2	Attach all LLD test reports: - Rack area - Dilution area		
9.3	Attach all QC kit results: 8x100 100&50 - Dilution - Val-Wash		
10	Accessories Manuals		
10.1	Check the manuals are available for the following accessories: - PC - Monitor - Printer - UPS (if present)		
11	Instrument Manuals		
11.1	Check the user manual availability and track: - Revision index - Language		
11.2	Check the troubleshooting guide availability and track: - Revision index - Language		
12	Comment / Observation:		

Engineer Name: _____ Signature: _____ Date: _____