

GRIFOLS	Reference	CTS-TD-SVD-000001	Version	8.0
	Status	Effective	Effective Date	23-Jul-2025
Printing Date	08-Aug-2025 06:18 UTC		Page	2 of 11
Title	Installation and Operational Qualification Report of the Procleix Panther System			

Purpose

This document provides a report to manually document the Installation and Operational Qualification of the Procleix Panther System or Procleix Panther System featuring ART. Alternatively, trained, qualified service personnel may opt to use the report generated through the Panther Service Manual to document IOQ activities.

The procedure follows the use of the Installation and Operational Qualification (IOQ) detailed in the Panther Service Manual. The use of this document shall be limited to trained, qualified service personnel when installing a Procleix Panther System or Procleix Panther System featuring ART. Service personnel performing Installation and Operational Qualification shall be trained on all the procedures as described in the Panther Service Manual.

References

- Procleix Panther System Service Manual
- Procleix Panther System Operator's Manual
- CTS-SRV-SVD-000032CTS-SRV-SVD-000032, Procleix Panther System Pre-site Evaluation Packet (completed copy for site)
- CTS-SRV-SVD-000033CTS-SRV-SVD-000033, Procleix Panther System Performance Qualification

Materials Required

Procleix Panther System or Procleix Panther System featuring ART

- Panther Instrument System
- Current Panther System Software, Firmware, Service Software, Teach Software
- Current Procleix Panther Assay Software
- Current Panther System Shield
- Current Panther System, Workstation Image, Image Deployment Tool

Laptop computer (Service Laptop)

- Current Panther Service Manual

GRIFOLS	Reference	CTS-TD-SVD-000001	Version	8.0
	Status	Effective	Effective Date	23-Jul-2025
Printing Date	08-Aug-2025 06:18 UTC		Page	3 of 11
Title	Installation and Operational Qualification Report of the Procleix Panther System			

Procedure

Perform the tasks associated with the appropriate Installation and Operational Qualification sections below as dictated by the configuration being installed. This document is to be used to document the Qualification activities as defined in the current version of the Panther Service Manual and DCM CTS-SRV-SVD-000033, Procleix Panther System Performance Qualification. The Procleix Panther System Operator's Manual can also be referenced for required operator specific tasks.

Complete the Installation and Qualification Checklist below, provide signature and date, and secure the customer signature and date on the Report. Include all created and associated data and documentation (files / printouts) with the final report in the CRM (Customer Relationship Management) application database Installation Record [Case or Work Order (WO)].

Documentation

All observations, results, deviations, and corrective actions will be documented in the Installation and Operational Qualification Report below. Reference the Procleix Panther System Operator's Manual and Service Manual as needed. The completed document shall provided to the customer. A copy is attached to the CRM Record or forwarded to regional Grifols Technical Service offices for attachment to the CRM Record.

GRIFOLS	Reference	CTS-TD-SVD-000001	Version	8.0
	Status	Effective	Effective Date	23-Jul-2025
Printing Date	08-Aug-2025 06:18 UTC		Page	4 of 11
Title	Installation and Operational Qualification Report of the Procleix Panther System			

Panther Installation and Operational Qualification Checklist

Site Name:				Date:	
Work Order #:		FSE:		Panther Serial #:	

Prerequisites

Review Pre-Site Evaluation Packet:

- All Pre-Site Evaluation Packet steps have been completed
- All Supplies have been received
- Any open items have been addressed

Parts and Materials Required:

- For high altitude sites (over 2000 feet or 610 meters), make sure that you have two new vacuum-sealed Auto Detect LCM Pumps available in case of persistent OLV errors.
- FSE Tool Kit
- Input Queue Teach Tool
- Shock Compression Tool for front panel springs
- TCR Carousel Teach Tool
- Reagent Bay, Pipettor Teach Tool (2 units needed for Full Auto Teach)
- Pipettor Teach Caps
- Cooling Module Fill Tube

For MTU Expansion

- MTU Queue Teach tool

All procedures can be found in the Panther System Service Manual **[MAN-06020]**, unless specified.

1 Uncrate & Inspect the Panther Instrument

- ☐ Uncrate the Instrument, the Accessory Kit, and the Marketing Panel.
- ☐ Confirm all Shock Watch/Tip-n-Tell Sensors are intact.
- ☐ Visually inspect and confirm ALL the Instrument's internal components are intact.

2 Instrument Placement

- ☐ Raise the feet of the instrument before moving it.

Moving the instrument with the feet lowered can bend the stabilizing rods.

- ☐ Move and place the Panther at the customer's specified location as per the Pre-Site Evaluation Packet.
- ☐ Set Up the Uninterrupted Power Supply (UPS). (Optional)

GRIFOLS	Reference	CTS-TD-SVD-000001	Version	8.0
	Status	Effective	Effective Date	23-Jul-2025
Printing Date	08-Aug-2025 06:18 UTC		Page	5 of 11
Title	Installation and Operational Qualification Report of the Procleix Panther System			

- ☐ Set the Input Voltage Range for the UPS. (Optional)
- ☐ Install the Monitor Arm and Monitor.
- ☐ Perform PC hardware setup: Connect all cables. (Follow section "Dell OptiPlex XE3 Installation" in Service Manual)

3 Software Installation

- ☐ Image the Panther PC.
- ☐ Reset the PC Workstation Clock and Date and configure as per customer preference (24- or 12-hour time).
- ☐ Change the PC Name.
- ☐ Install Software as per relevant FSB and Service Manual.
- ☐ Generate Installed Assays Report as per relevant FSB and attach to WO.

4 Initial Module Preparation/Instrument Preparation

- ☐ Remove all packing foam on the Pipettor Gantry.
- ☐ Open the Service Drawer and remove all packing foam blocks from around the Linear Distributor.
- ☐ Check to be sure that all modules are intact and have not shifted during relocation.
- ☐ Check all modules for missing screws and that all screws are tight. Install the Side Panels, Canopies and Front Canopy Doors.
- ☐ Install the Handheld Barcode Scanner Bracket.
- ☐ Verify the MTU Input Queue ball detent is correctly adjusted and lubricated [SAP 741459/MME-02134 Grease, Linear Bearings] *If the ball detent is not lubed, apply grease.*
- ☐ Inspect and verify that all three incubator lids are fully seated.
- ☐ Inspect and verify that all 6 rubber restraints are on all incubators and are not damaged.
- ☐ Check incubator doors for sticking – *lightly* lube door rail with [SAP717920/903864 Standard Lubrication Kit] if door is sticky.
- ☐ Install Peristaltic Pump tubing. *Ensure the Luer fittings are tight. Do NOT over torque the fittings. Use tape and label the pump with the date of installation.*
- ☐ Open the Waste Drawer and apply grease [SAP716265/903524 Dow Corning High Vacuum Grease] to the Waste Bottle fittings and O-rings.
- ☐ Retrieve the Output Queue Lid from the Universal Fluids Drawer and install it on the Output Queue.
- ☐ Close the Service Drawer and route tubing and wires so that they do not interfere with the Linear Distributor or Pipettor movement. Check the left back corner where the Service Drawer rail touches the back of the Panther.
- ☐ Hand-tighten the two Pipettor DiTi thumbscrews.
- ☐ Level the instrument. Ensure that the TCR Door closes without additional friction and that it is aligned.

GRIFOLS	Reference	CTS-TD-SVD-000001	Version	8.0
	Status	Effective	Effective Date	23-Jul-2025
Printing Date	08-Aug-2025 06:18 UTC		Page	6 of 11
Title	Installation and Operational Qualification Report of the Procleix Panther System			

4B Use Blue sections only for installations including MTU Expansion

- ☐ Install new MTU sidecar CAN Y-splitter cable
- ☐ Remove Input Queue Install MTU Queue
- ☐ Move MTU Side car close to the Panther and connect CAN splitter cable
- ☐ Run Instrument Setup – select all relevant options including Panther, Waste on the Fly and MTU Expansion
- ☐ Verify Distributor Theta Alignment.
- ☐ Align the Distributor to the MTU Queue.
- ☐ Align the Distributor to the MTU Scanner.
- ☐ Using MTU Teach Tool, Align and Teach Linear Distributor to MTU Queue
- ☐ Install guide plate
- ☐ Rough Align MTU Sidecar with power off
- ☐ Align Linear Distributor to MTU Queue
- ☐ Set guide plate
- ☐ Align MTU Sidecar.
- ☐ Perform Full MTU Expansion Autoteach (Lift)
- ☐ Perform End of Travel teaching for the pusher
- ☐ Install Carriers
- ☐ Install the Handheld Barcode Scanner Bracket.

4C Use Red sections only for installations including Waste to Drain

- ☐ Inspect customer’s drain for obstructions or clogging
- ☐ If Panther is Serial #1936 and below, punch or drill hole in rear of Panther
- ☐ Drill hole in Waste Drawer for waste to drain tubing
- ☐ Remove plug from waste bottle and install
- ☐ Install Pump Module Install PTD Tubing
- ☐ Install Power and Sensor Cables and Grommet

GRIFOLS	Reference	CTS-TD-SVD-000001	Version	8.0
	Status	Effective	Effective Date	23-Jul-2025
Printing Date	08-Aug-2025 06:18 UTC		Page	7 of 11
Title	Installation and Operational Qualification Report of the Procleix Panther System			

5 Run Instrument Setup (Firmware)

☐ Run Instrument Setup (Loading Firmware). Make sure to choose the correct options on the Setup Splash Screen.

Never select Fusion, RTFs, or Sidecar.

These are the options:

- APP Pipettors or BLDC Pipettors (select the type you are installing)
- MTU Expansion (yes/no)
- Continuous Access Waste and Universal Fluids (yes/no)
 - o Pump to Drain (Waste to Drain) (yes/no)
- Track Connectivity (yes/no)

6 Verify ALL fittings are hand-tight

- ☐ Oil, Wash, Auto Detect 1, and Auto Detect 2 pump fittings.
- ☐ Universal Fluid Drawer solenoid valve fittings.
- ☐ MagWash 1 and MagWash 2 valve manifold fittings.
- ☐ Magwash Aspirators (use a 7mm wrench).
- ☐ Peri-Pump fittings.
- ☐ Vacuum Manifold fittings.
- ☐ Vacuum Pressure Sensor fittings.
- ☐ Oil Pressure Sensor fittings.
- ☐ Universal Fluids Drawer Oil and Wash 3-way manifold fittings.
- ☐ Make sure all drawers are closed completely and securely.
- ☐ Install the Waste Bottle, Waste Bags, and Waste Deflector
- ☐ Install the Handheld Barcode Scanner.

7 Pipettor Auto-Teaching

- ☐ Make a backup of the coordinate.coo file before starting the Auto-Teaching.
- ☐ Teach ALL Modules (Pipettor Teaching, Automatic (All Modules)).

GRIFOLS	Reference	CTS-TD-SVD-000001	Version	8.0
	Status	Effective	Effective Date	23-Jul-2025
Printing Date	08-Aug-2025 06:18 UTC		Page	8 of 11
Title	Installation and Operational Qualification Report of the Procleix Panther System			

8 Linear Distributor Alignment

- ☐ Verify free movement in all Distributor axis.
- ☐ Align Theta position in Linear Distributor and Lumo.
- ☐ Perform MTU Input Queue Adjustment in the Y-Axis.
- ☐ Align the Distributor to the MTU Input Queue with or without Input Queue Teaching Tool.
- ☐ Align the Distributor to the MTU Scanner.
- ☐ Execute Full Auto-Teach.
- ☐ Perform manual Distributor Alignment to the Service Queue.

9 Pipettor Pressure Integrity Test

- ☐ Perform the Pipettor Integrity Pressure Test on the Reagent Pipettor Arm.
- ☐ Perform the Pipettor Integrity Pressure Test on the Sample Pipettor Arm.

10A Service Software Checks and OQs

- ☐ Fill the Cooling Module.
- ☐ Configure the Reagent and Sample Bay Barcode Scanners.
- ☐ Setup the Service Software Report Header.
- ☐ Perform 2 cycles of System Level OQ Procedure, selecting all modules.

10B Use Green sections only for installations including CONTINUOUS ACCESS WASTE & FLUIDS

- ☐ Perform 2 cycles of Solid Waste on the Fly Verification Procedure/OQ in Service Software
- ☐ Perform Liquid Waste on the fly Verification Procedure

11 Final Checks

- ☐ Install the Pipettor Front Shield.
- ☐ Install the Front Marketing Panel.
- ☐ Install the Trims and Door Panels.

GRIFOLS	Reference	CTS-TD-SVD-000001	Version	8.0
	Status	Effective	Effective Date	23-Jul-2025
Printing Date	08-Aug-2025 06:18 UTC		Page	9 of 11
Title	Installation and Operational Qualification Report of the Procleix Panther System			

12 System Bootup

- ☐ Load MTUs into the MTU Input Queue.
- ☐ Load Tips into the Tip Drawers.
- ☐ Install Fluids AD1, AD2, Oil, MagWash Buffer, Bleach, and Deactivation Solution.
- ☐ Install the Universal Fluid Bottle connectors.
- ☐ Open the Panther Main software and verify that the instrument boots without errors.
- ☐ Update the Instrument Serial Number in Operator Software.
- ☐ Create an Administrator Account in Operator Software.
- ☐ Select the System Software Language.
- ☐ Install the appropriate model Printer or Network Printer according with service manual.

If not a supported printer, search for brand and model during installation procedure, if model is not available choose Microsoft PLC6 Class Driver for USB Printers, verify printer installs and prints a test page and a report in Panther GUI.

- ☐ Shut down and Restart the Instrument and PC. **Allow the PC to shut down completely. Do not force a shutdown.**
- ☐ Complete two Full Primes.
- ☐ Inspect ALL fluid lines for leaks.
- ☐ Inspect Universal Fluid Drawer for leaks.

13 Setup the Firewall and Network Connections (if needed)

- ☐ Configure the firewall with the latest version of Firewall Wizard.
- ☐ Connect the configured firewall to the Panther Instrument.
- ☐ Configure the Panther Workstation with the appropriate static IP address.
- ☐ If applicable, install and enable Secure Remote Connectivity software (SRC) as per CTS-TD-SVD-000008 Grifols Secure Remote Connectivity Installation and Usage Field Support (Optional)
- ☐ If applicable, configure Laboratory Information System (LIS) connection. Select **ONE**:
 - ☐ LIS Serial Install **OR** ☐ LIS TCP/IP Install

14 Panther Luminometer Flashcheck Reagent Verification Procedure

- ☐ Perform the Flashcheck Reagent Verification procedure.
- ☐ Save the report to a flash drive for attachment to the Installation work order.

GRIFOLS	Reference	CTS-TD-SVD-000001	Version	8.0
	Status	Effective	Effective Date	23-Jul-2025
Printing Date	08-Aug-2025 06:18 UTC		Page	10 of 11
Title	Installation and Operational Qualification Report of the Procleix Panther System			

15 Panther System Performance Qualification (PQ)

- ☐ Perform required Performance Qualification(s) as per CTS-SRV-SVD-000033 Procleix Panther System Performance Qualification
- ☐ Save the report to a flash drive for attachment to the Installation Work order.

16 Post Qualification

- ☐ Dispose of any partial assay reagents according to laboratory procedures.
- ☐ Empty solid and liquid waste according to laboratory procedures

17 Additional documentation.

- ☐ Copy of the Pre-Site Evaluation (To be discussed if it's useful/necessary)
- ☐ The System Level OQ Report
- ☐ The Pipettor Integrity Pressure test
- ☐ The FlashCheck results page
- ☐ The Performance Qualification results
- ☐ Installed Assays Report verified against latest relevant Assay Installer. Identify the "Intended use", CE-IVD or US-IVD/IQC, writing it down on the printed report. Add a check mark in each assay version from the report once verified against the corresponding one in the Assay Installer table of the relevant Technical/ Field Service Bulletin. Date and sign the report.

