

Document Detail

12/19/2017

1

Type: Form
 Document No.: FRM-471908650F1[B]
 Title: AmiCORE Installation and Preventative Maintenance Checklist
 Owner: 180305 Jeffery Mashke
 Status: CURRENT
 Effective Date: 19-Dec-2017 Expiration Date: 19-Dec-2019

Reference

<u>Document No.</u>	<u>Content Type</u>	<u>Relation</u>	<u>Fixed Rev</u>	<u>Status</u>
SDL-SDL12830001[K]	DOCUMENT	Master Document	No	CURRENT
Title:AmiCORE Service Document Listing (SDL)				

Activities

<u>Document Build No.</u>	<u>Access Activity</u>	<u>Accessed By</u>	<u>Accessed Date</u>
1	Check In	180305	13-Dec-2017

Note:

Review

Build No.: 1

Closed Date: 12/19/2017 4:58:29PM

Review: Document Approval

Review Purpose:

Review Note: SYSTEM AUTO CLOSE REVIEW

<u>Level</u>	<u>Owner Role</u>	<u>Actor</u>	<u>Sign-off By</u>
1	Document Owner Document Owner	180305 Jeffery Mashke	180305

Note To Approver:

Note From Approver:approved

2	Document Control Document Control	979631 Susan Scheef	979631
---	-----------------------------------	---------------------	--------

Note To Approver:

Note From Approver:approved

1	Quality Approver Quality Approver	987007 Peter Kim	987007
---	-----------------------------------	------------------	--------

Note To Approver:

Note From Approver:Reviewed and approved

AmiCORE Installation/Preventive Maintenance Checklist

AmiCORE S/N: _____
 Customer: _____
 Service #: _____

Date: _____
 SW Version: _____
 FSE Name: _____

Instructions:

1. For detailed information, including required test points, refer to the AmiCORE Service Manual.
2. For installations, complete Sections I and III. For Preventative Maintenance, complete Sections II and III.

Check box note: P = Pass, F = Fail, N/A = Not Applicable

I. Installation Only		NA ☐
3.2	Instrument Unpacking & Assembly	P ☐ F ☐
6.1	Touch Screen Calibration	P ☐ F ☐
6.2	Weight Scale Verification & Calibration	P ☐ F ☐
6.3	Centrifuge Pressure Sensor Calibration	P ☐ F ☐
6.4	Interface Detector Assembly Alignment Calibration	P ☐ F ☐
II. Preventative Maintenance Only		NA ☐
4.2	Parts Replacement P.M.	P ☐ F ☐
III. Installation & Preventive Maintenance		
3.3	Visual Inspection	Completed ☐
3.4	Power UP System Test	P ☐ F ☐
3.5	Enable Service Mode	P ☐ F ☐
3.6	Pump Verification	P ☐ F ☐
3.7	Cassette Motor Verification	P ☐ F ☐
3.8	Clamp Verification	P ☐ F ☐
3.9	Valve Verification	P ☐ F ☐
3.10	Pressure Sensor Verification	P ☐ F ☐
3.11	Air Detector Verification	P ☐ F ☐
3.12	Door Lock Verification	P ☐ F ☐
3.13	Centrifuge Verification	P ☐ F ☐
3.14	Interface Verification	P ☐ F ☐
3.15	Weight Scale Verification	P ☐ F ☐
3.16	Leak Detector Verification	P ☐ F ☐
3.17	Temperature Sensor Verification	P ☐ F ☐
3.18	Cuff Verification	P ☐ F ☐
3.19	USB Verification	P ☐ F ☐
3.20	Power Verification	P ☐ F ☐
3.21	Battery Verification	P ☐ F ☐
3.22	Simulated AmiCORE Procedure	P ☐ F ☐

Comments: _____

Signature: _____ Date: _____