

LabUMat Installation Checklist

Work Order No.: Installation Date:

Customer Name: Department:

Instrument Serial No..... Instrument SW Version :.....

No	Procedure	Status
1	Unpack the Instrument and place it on the table.	
2	Remove the pipetting probe retainer plate.	
3	Install the liquid tanks.	
4	Install the rack conveyor, comb conveyor and pipetting stage.	
5	Check and adjust if necessary the position of the strip clamps.	
6	Switch on the device and initialize.	
7	Check and adjust the front-back axis alignment of the linear robot horizontal motor.	
8	Check and adjust if necessary the pipetting position.	
9	Install The rack conveyor and check the rack motion is smooth.	
10	Check and adjust if necessary the pad-reading position.	
11	Check and adjust if necessary the rack catch corrections.	
12	Set up the Flipped strip sensor.	
13	Set up the test tube minimum sample level.	
14	Load strips and run sample/control to check overall functioning	

Installed By: Er. Signature: