

LabUMat Maintenance Checklist

Cutomer :	S\N:				
Department :	Average throughput of the instrument (measurement / day):				
Date :	1-150 meas./day ☐ 150-300 meas./day ☐ 300+ meas./day ☐				
Feeder module	Time req'd	1-150 meas./day	150-300 meas./day	300+ meas./day	done
Clean the surface of the feeder cylinder and the strip sliding ramp	3 min.	every 6 months	every 4 months	every 2 month	☐
Grease the side of the feeder cylinder	2 min.				☐
Check the cleanliness and the reflectance values of the strip presence&orientation sensor	3 min.				☐
Check the position of feeder blade and readjust if necessary	5 min.	every 6 months	every 4 months	every 2 months	☐
Check the pipetting position of strip pusher and readjust if necessary	3 min.				☐
Check the operation of all sensors and motors	5 min.	every 12 months	every 12 months	every 6 months	☐
Wash module	Time req'd	1-150 meas./day	150-300 meas./day	300+ meas./day	done
Perform a dry test	5 min.	every 6 months	every 4 months	every 2 months	☐
Replace the components included in the maintenance kit	15 min.	every 12 months	every 12 months	every 6 months	☐
Robot module	Time req'd	1-150 meas./day	150-300 meas./day	300+ meas./day	done
Check the straightness of the pipette and the smooth lowering into the washing house	3 min.	every 6 months	every 4 months	every 2 months	☐
Check the pipetting positions	3 min.				☐
Check the operation of all sensors and motors	5 min.	every 12 months	every 12 months	every 12 months	☐
Rack in-out module	Time req'd	1-150 meas./day	150-300 meas./day	300+ meas./day	done
Check rack mover height	2 min.	every 6 months	every 4 months	every 2 months	☐
Check the tongue parameters	5 min.				☐
Clean the gear belt of the rack in/out unit	5 min.				☐
Clean the rack mover inside and outside	8 min.				☐
Clean and lubricate the gear belts of the rack mover unit	3 min.				☐
Check the operation of all sensors and motors	5 min.	every 12 months	every 12 months	every 6 months	☐
Strip-Head moving module	Time req'd	1-150 meas./day	150-300 meas./day	300+ meas./day	done
Check the Head Strip position	3 min.	every 6 months	every 4 months	every 2 months	☐
Check the Reference Pad position	3 min.				☐
Replace the components included in the maintenance kit	10 min.	every 12 months	every 12 months	every 6 months	☐
Check the operation of all sensors and motors	5 min.				☐
Sampler module	Time req'd	1-150 meas./day	150-300 meas./day	300+ meas./day	done
Perform a wet test for the aspiration system	3 min.	every 6 months	every 4 months	every 2 months	☐
Replace the components included in the maintenance kit	2 min.	every 12 months	every 12 months	every 6 months	☐
Check the operation of all sensors and motors	5 min.				☐
Power supply	Time req'd	1-150 meas./day	150-300 meas./day	300+ meas./day	done
Clean the inside with pressurized air	2 min.	every 12 months	every 12 months	every 6 months	☐
Embedded PC	Time req'd	1-150 meas./day	150-300 meas./day	300+ meas./day	done
Perform Disk Cleanup, Error Checking and Defragmentation	10 min.	every 12 months	every 12 months	every 6 months	☐

Note: Required time does not include the time for removing the back panel and removing the module or part (if necessary)

Service Engineer :	Customer Signatuer :
---------------------------	-----------------------------