



UriSed_Installation Check list

Work Order No.: Installation Date:

Customer Name: Department:

Instrument Serial No..... Instrument SW Version :.....

No	Procedure	Status
1	Unpack the Instrument and place it on the table.	
2	Unscrew the retainer bolt that secures the robot during transport.	
3	Place the rotor on the pin of the axis holder.	
4	Screw the rotor fixing and put the feeder top disk.	
5	Insert the front cuvette rail & back cuvette rail.	
6	Install the fluidic system.	
7	Install the rack conveyor.	
8	Install the PC and its peripherals.	
9	Switch on the device and Initialize.	
10	If the device is connected with LabUMat check the connectivity.	
11	Check and adjust if necessary rack puller correction.	
12	Check and adjust if necessary tube puller correction.	
13	Check and adjust if necessary the stub position.	
14	Check and adjust if necessary the cuvette detection laser position.	
15	Check and adjust if necessary the centrifuge null position.	
16	Review the microscope set up.	
17	Check and adjust if necessary the pipetting position.	
18	Check and adjust the rack catch correction.	
19	Setup the test tube minimum sample level.	
20	Run sample/control to check overall functioning	
21	run empty cuvette verify and run QC	

Installed By: Er. Signature: