

UriSed Maintenance Checklist

Hospital :	S\N:				
Department :	Average throughput of the instrument (measurement / day):				
Date:	1-150 meas./day ☐ 150-300 meas./day ☐ 300+ meas./day ☐				
Feeder module	Time req'd	1-150 meas./day	150-300 meas./day	300+ meas./day	done
Check the pushing arm + cuvette contact point	3 min.	every 6 months	every 4 months	every 2 month	☐
Check the Pipetting position of the arm	5 min.				☐
Check the cleanliness of the front cuv.rail pipe	3 min.	every 6 months	every 4 months	every 2 months	☐
Check the snap-fits of the arm	3 min.				☐
Check the centrifuge and microscope pos. of the arm	5 min.	every 12 months	every 12 months	every 6 months	☐
Check the condition of the cuvette front rail spring	5 min.				☐
Check the operation of all sensors and motors	5 min.				☐
Wash module	Time req'd	1-150 meas./day	150-300 meas./day	300+ meas./day	done
Perform a dry test	3 min.	every 6 months	every 4 months	every 2 months	☐
Replace the components included in the maintenance kit	15 min.	every 12 months	every 12 months	every 6 months	☐
Robot module	Time req'd	1-150 meas./day	150-300 meas./day	300+ meas./day	done
Check the condition of the pipette (straightness, coating)	3 min.	every 6 months	every 4 months	every 2 months	☐
Check the pipetting position and filling process	3 min.				☐
Perform a wet test for the aspiration system	3 min.				☐
Replace the components included in the maintenance kit	10 min.	every 12 months	every 12 months	every 6 months	☐
Check the operation of all sensors and motors	5 min.				☐
Rack in-out module	Time req'd	1-150 meas./day	150-300 meas./day	300+ meas./day	done
Check rack mover height	2 min.	every 6 months	every 4 months	every 2 month	☐
Check the tongue parameters	5 min.				☐
Clean the gear belt of the rack in/out unit	5 min.	every 6 months	every 4 months	every 2 months	☐
Clean the rack mover inside and outside	8 min.				☐
Clean and lubricate the gear belts of the rack mover unit	3 min.				☐
Check the operation of all sensors and motors	5 min.	every 12 months	every 12 months	every 6 months	☐
Microscope module	Time req'd	1-150 meas./day	150-300 meas./day	300+ meas./day	done
Check the cleanliness of the objective lens and the lamp unit	1 min.	every 6 months	every 4 months	every 2 months	☐
Check the focus graph with the reference cuvette	10 min.				☐
Check the operation of all sensors and motors	5 min.	every 12 months	every 12 months	every 6 months	☐
Centrifuge module	Time req'd	1-150 meas./day	150-300 meas./day	300+ meas./day	done
Check the zero position of the centrifuge arm	3 min.	every 6 months	every 4 months	every 2 months	☐
Check the cleanliness of the inside of the centrifuge cover	2 min.	every 12 months	every 12 months	every 6 months	☐
Power supply	Time req'd	1-150 meas./day	150-300 meas./day	300+ meas./day	done
Clean the inside with pressurized air	2 min.	every 12 months	every 12 months	every 6 months	☐
PC	Time req'd	1-150 meas./day	150-300 meas./day	300+ meas./day	done
Perform Disk Cleanup, Error Checking and Defragmentation	10 min.	every 12 months	every 12 months	every 6 months	☐

Note: Required time does not include the time for removing the back panel and removing the module or part (if necessary)

Service Engineer :	Customer Signatuer :
--------------------	----------------------