

RESULTS BY WORKLIST

Serial Number:	00126	System SW Version:	3.2.1.6	Date Range:	3/25/2011 15:22:10 to 3/25/2011 19:17:09
Assay:	CT/GC				
Worklist IDs:	00126-20110325-03 Version: 1.33 (For In Vitro Diagnostic Use) ML: 585289 Exp: 6/15/2012				
Operators:	fhaider				
Specimen Tested:	173	GC POS: 112	CT POS: 4	GC EQUIV: 3	CT EQUIV: 1
Status Flag Legend:	M - Sample ID of a new test request was manually entered via the keyboard when a barcode was not readable by the instrument.				
Reviewed By:		Date:		Verified By:	Date:

CT/GC

Version: 1.33 (For In Vitro Diagnostic Use)

Worklist ID: 00126-20110325-03

Time	Rack-Tube #	Sample ID	Total RLU	CT Result	GC Result	Flags
3/25/2011 15:22:10	835-02	CT Positive Control - Lot 586950	838	CT POS	GC neg	
3/25/2011 15:22:10	835-01	GC Positive Control - Lot 586951	1124	CT neg	GC POS	
3/25/2011 15:22:11	835-03	0000001	16	CT neg	GC neg	M
3/25/2011 15:22:11	835-04	0000002	6	CT neg	GC neg	M
3/25/2011 15:22:11	835-05	0000003	59	CT neg	GC neg	M
3/25/2011 15:22:11	835-06	0000004	5	CT neg	GC neg	M
3/25/2011 15:27:10	835-07	0000005	4	CT neg	GC neg	M
3/25/2011 15:27:10	835-08	0000006	4	CT neg	GC neg	M
3/25/2011 15:27:10	835-09	0000007	2	CT neg	GC neg	M
3/25/2011 15:27:11	835-10	0000008	4	CT neg	GC neg	M
3/25/2011 15:32:10	835-11	0000009	5	CT neg	GC neg	M
3/25/2011 15:32:10	835-12	0000010	4	CT neg	GC neg	M
3/25/2011 15:32:10	835-13	0000011	4	CT neg	GC neg	M
3/25/2011 15:32:10	835-14	0000012	4	CT neg	GC neg	M

CT/GC

Version: 1.33 (For In Vitro Diagnostic Use)

Time	Rack-Tube #	Sample ID	Total RLU	CT Result	GC Result	Flags
3/25/2011 15:32:10	835-15	0000013	4	CT neg	GC neg	M
3/25/2011 15:37:10	945-01	0000014	5	CT neg	GC neg	M
3/25/2011 15:37:10	945-02	0000015	5	CT neg	GC neg	M
3/25/2011 15:37:10	945-03	0000016	4	CT neg	GC neg	M
3/25/2011 15:37:10	945-04	0000017	4	CT neg	GC neg	M
3/25/2011 15:37:10	945-05	0000018	4	CT neg	GC neg	M
3/25/2011 15:42:12	945-06	0000019	5	CT neg	GC neg	M
3/25/2011 15:42:12	945-07	0000020	4	CT neg	GC neg	M
3/25/2011 15:42:12	945-08	0000021	5	CT neg	GC neg	M
3/25/2011 15:42:12	945-09	0000022	3	CT neg	GC neg	M
3/25/2011 15:42:12	945-10	0000023	4	CT neg	GC neg	M
3/25/2011 15:47:10	945-11	0000024	5	CT neg	GC neg	M
3/25/2011 15:47:10	945-12	0000025	5	CT neg	GC neg	M
3/25/2011 15:47:10	945-13	0000026	4	CT neg	GC neg	M
3/25/2011 15:47:10	945-14	0000027	3	CT neg	GC neg	M
3/25/2011 15:47:10	945-15	0000028	4	CT neg	GC neg	M
3/25/2011 15:52:12	848-01	0000029	5	CT neg	GC neg	M
3/25/2011 15:52:12	848-02	0000030	5	CT neg	GC neg	M
3/25/2011 15:52:12	848-03	0000031	4	CT neg	GC neg	M
3/25/2011 15:52:13	848-04	0000032	5	CT neg	GC neg	M
3/25/2011 15:52:13	848-05	0000033	4	CT neg	GC neg	M
3/25/2011 15:57:12	848-06	0000034	5	CT neg	GC neg	M
3/25/2011 15:57:12	848-07	0000035	4	CT neg	GC neg	M
3/25/2011 15:57:12	848-08	0000036	4	CT neg	GC neg	M
3/25/2011 15:57:13	848-09	0000037	4	CT neg	GC neg	M
3/25/2011 15:57:13	848-10	0000038	3	CT neg	GC neg	M
3/25/2011 16:02:10	848-11	0000039	5	CT neg	GC neg	M
3/25/2011 16:02:10	848-12	0000040	5	CT neg	GC neg	M
3/25/2011 16:02:10	848-13	0000041	4	CT neg	GC neg	M
3/25/2011 16:02:10	848-14	0000042	4	CT neg	GC neg	M
3/25/2011 16:02:10	848-15	0000043	3	CT neg	GC neg	M
3/25/2011 16:07:09	947-01	0000044	5	CT neg	GC neg	M

CT/GC

Version: 1.33 (For In Vitro Diagnostic Use)

Time	Rack-Tube #	Sample ID	Total RLU	CT Result	GC Result	Flags
3/25/2011 16:07:10	947-02	0000045	5	CT neg	GC neg	M
3/25/2011 16:07:10	947-03	0000046	4	CT neg	GC neg	M
3/25/2011 16:07:10	947-04	0000047	4	CT neg	GC neg	M
3/25/2011 16:07:10	947-05	0000048	4	CT neg	GC neg	M
3/25/2011 16:12:10	947-06	0000049	5	CT neg	GC neg	M
3/25/2011 16:12:10	947-07	0000050	5	CT neg	GC neg	M
3/25/2011 16:12:10	947-08	0000051	4	CT neg	GC neg	M
3/25/2011 16:12:10	947-09	0000052	5	CT neg	GC neg	M
3/25/2011 16:12:10	947-10	0000053	1177	CT neg	GC POS	M
3/25/2011 16:17:10	947-11	0000054	1104	CT neg	GC POS	M
3/25/2011 16:17:10	947-12	0000055	1054	CT neg	GC POS	M
3/25/2011 16:17:10	947-13	0000056	1046	CT neg	GC POS	M
3/25/2011 16:17:10	947-14	0000057	1115	CT neg	GC POS	M
3/25/2011 16:17:10	947-15	0000058	1121	CT neg	GC POS	M
3/25/2011 16:22:10	849-01	0000059	1189	CT neg	GC POS	M
3/25/2011 16:22:10	849-02	0000060	1112	CT neg	GC POS	M
3/25/2011 16:22:10	849-03	0000061	1125	CT neg	GC POS	M
3/25/2011 16:22:10	849-04	0000062	2078	CT POS	GC POS	M
3/25/2011 16:22:10	849-05	0000063	643	CT POS	GC POS	M
3/25/2011 16:27:10	849-06	0000064	680	CT neg	GC POS	M
3/25/2011 16:27:10	849-07	0000065	992	CT neg	GC POS	M
3/25/2011 16:27:10	849-08	0000066	977	CT neg	GC POS	M
3/25/2011 16:27:10	849-09	0000067	1271	CT POS	GC POS	M
3/25/2011 16:27:10	849-10	0000068	811	CT neg	GC POS	M
3/25/2011 16:32:10	849-11	0000069	774	CT neg	GC POS	M
3/25/2011 16:32:10	849-12	0000070	254	CT neg	GC POS	M
3/25/2011 16:32:10	849-13	0000071	25	CT neg	GC neg	M
3/25/2011 16:32:10	849-14	0000072	416	CT neg	GC POS	M
3/25/2011 16:32:11	849-15	0000073	484	CT neg	GC POS	M
3/25/2011 16:37:09	839-01	0000074	250	CT neg	GC POS	M
3/25/2011 16:37:10	839-02	0000075	827	CT neg	GC POS	M
3/25/2011 16:37:10	839-03	0000076	125	CT neg	GC EQUIV	M

CT/GC

Version: 1.33 (For In Vitro Diagnostic Use)

Time	Rack-Tube #	Sample ID	Total RLU	CT Result	GC Result	Flags
3/25/2011 16:37:10	839-04	0000077	180	CT neg	GC POS	M
3/25/2011 16:37:10	839-05	0000078	309	CT neg	GC POS	M
3/25/2011 16:42:09	839-06	0000079	337	CT neg	GC POS	M
3/25/2011 16:42:10	839-07	0000080	215	CT neg	GC POS	M
3/25/2011 16:42:10	839-08	0000081	270	CT neg	GC POS	M
3/25/2011 16:42:10	839-09	0000082	48	CT neg	GC neg	M
3/25/2011 16:42:10	839-10	0000083	210	CT neg	GC POS	M
3/25/2011 16:47:10	839-11	0000084	154	CT neg	GC POS	M
3/25/2011 16:47:10	839-12	0000085	24	CT neg	GC neg	M
3/25/2011 16:47:10	839-13	0000086	10	CT neg	GC neg	M
3/25/2011 16:47:10	839-14	0000087	15	CT neg	GC neg	M
3/25/2011 16:47:10	839-15	0000088	26	CT neg	GC neg	M
3/25/2011 16:52:09	946-01	0000089	891	CT neg	GC POS	M
3/25/2011 16:52:10	946-02	0000090	1072	CT neg	GC POS	M
3/25/2011 16:52:10	946-03	0000091	1141	CT neg	GC POS	M
3/25/2011 16:52:10	946-04	0000092	737	CT neg	GC POS	M
3/25/2011 16:52:10	946-05	0000093	621	CT neg	GC POS	M
3/25/2011 16:57:10	946-06	0000094	427	CT neg	GC POS	M
3/25/2011 16:57:10	946-07	0000095	880	CT neg	GC POS	M
3/25/2011 16:57:10	946-08	0000096	845	CT neg	GC POS	M
3/25/2011 16:57:10	946-09	0000097	2124	CT POS	GC POS	M
3/25/2011 16:57:10	946-10	0000098	866	CT neg	GC POS	M
3/25/2011 17:02:10	946-11	0000099	622	CT neg	GC POS	M
3/25/2011 17:02:10	946-12	0000100	592	CT neg	GC POS	M
3/25/2011 17:02:10	946-13	0000101	497	CT neg	GC POS	M
3/25/2011 17:02:10	946-14	0000102	458	CT neg	GC POS	M
3/25/2011 17:02:10	946-15	0000103	853	CT neg	GC POS	M
3/25/2011 17:07:09	836-01	0000104	1053	CT neg	GC POS	M
3/25/2011 17:07:10	836-02	0000105	217	CT neg	GC POS	M
3/25/2011 17:07:10	836-03	0000106	475	CT neg	GC POS	M
3/25/2011 17:07:10	836-04	0000107	919	CT neg	GC POS	M
3/25/2011 17:07:10	836-05	0000108	345	CT neg	GC POS	M

CT/GC

Version: 1.33 (For In Vitro Diagnostic Use)

Time	Rack-Tube #	Sample ID	Total RLU	CT Result	GC Result	Flags
3/25/2011 17:12:09	836-06	0000109	213	CT neg	GC POS	M
3/25/2011 17:12:10	836-07	0000110	206	CT neg	GC POS	M
3/25/2011 17:12:10	836-08	0000111	249	CT EQUIV	GC EQUIV	M
3/25/2011 17:12:10	836-09	0000112	189	CT neg	GC POS	M
3/25/2011 17:12:10	836-10	0000113	1245	CT neg	GC POS	M
3/25/2011 17:12:10	836-11	0000114	1221	CT neg	GC POS	M
3/25/2011 17:17:10	836-12	0000115	1100	CT neg	GC POS	M
3/25/2011 17:17:10	836-13	0000116	1141	CT neg	GC POS	M
3/25/2011 17:17:10	836-14	0000117	1170	CT neg	GC POS	M
3/25/2011 17:17:10	836-15	0000118	1106	CT neg	GC POS	M
3/25/2011 17:22:09	835-03	0000119	1147	CT neg	GC POS	M
3/25/2011 17:22:09	835-04	0000120	915	CT neg	GC POS	M
3/25/2011 17:22:09	835-05	0000121	1137	CT neg	GC POS	M
3/25/2011 17:22:09	835-06	0000122	1079	CT neg	GC POS	M
3/25/2011 17:22:10	835-07	0000123	1058	CT neg	GC POS	M
3/25/2011 17:27:09	835-08	0000124	1083	CT neg	GC POS	M
3/25/2011 17:27:09	835-09	0000125	1034	CT neg	GC POS	M
3/25/2011 17:27:09	835-10	0000126	1072	CT neg	GC POS	M
3/25/2011 17:27:09	835-11	0000127	1043	CT neg	GC POS	M
3/25/2011 17:27:09	835-12	0000128	955	CT neg	GC POS	M
3/25/2011 17:32:08	835-13	0000129	1033	CT neg	GC POS	M
3/25/2011 17:32:09	835-14	0000130	916	CT neg	GC POS	M
3/25/2011 17:32:09	835-15	0000131	1032	CT neg	GC POS	M
3/25/2011 17:32:09	945-01	0000132	864	CT neg	GC POS	M
3/25/2011 17:32:09	945-02	0000133	972	CT neg	GC POS	M
3/25/2011 17:37:08	945-03	0000134	964	CT neg	GC POS	M
3/25/2011 17:37:09	945-04	0000135	930	CT neg	GC POS	M
3/25/2011 17:37:09	945-05	0000136	802	CT neg	GC POS	M
3/25/2011 17:37:09	945-06	0000137	954	CT neg	GC POS	M
3/25/2011 17:37:09	945-07	0000138	868	CT neg	GC POS	M
3/25/2011 17:42:09	945-08	0000139	751	CT neg	GC POS	M
3/25/2011 17:42:09	945-09	0000140	750	CT neg	GC POS	M

CT/GC

Version: 1.33 (For In Vitro Diagnostic Use)

Time	Rack-Tube #	Sample ID	Total RLU	CT Result	GC Result	Flags
3/25/2011 17:42:10	945-10	0000141	864	CT neg	GC POS	M
3/25/2011 17:42:10	945-11	0000142	681	CT neg	GC POS	M
3/25/2011 17:42:10	945-12	0000143	900	CT neg	GC POS	M
3/25/2011 18:07:09	945-13	0000144	379	CT neg	GC POS	M
3/25/2011 18:07:10	945-14	0000145	411	CT neg	GC POS	M
3/25/2011 18:07:10	945-15	0000146	765	CT neg	GC POS	M
3/25/2011 18:07:10	848-01	0000147	541	CT neg	GC POS	M
3/25/2011 18:07:10	848-02	0000148	501	CT neg	GC POS	M
3/25/2011 18:32:09	848-03	0000149	864	CT neg	GC POS	M
3/25/2011 18:32:09	848-04	0000150	387	CT neg	GC POS	M
3/25/2011 18:32:09	848-05	0000151	131	CT neg	GC EQUIV	M
3/25/2011 18:32:09	848-06	0000152	706	CT neg	GC POS	M
3/25/2011 18:32:10	848-07	0000153	339	CT neg	GC POS	M
3/25/2011 18:37:09	848-08	0000154	392	CT neg	GC POS	M
3/25/2011 18:37:09	848-09	0000155	809	CT neg	GC POS	M
3/25/2011 18:37:10	848-10	0000156	844	CT neg	GC POS	M
3/25/2011 18:37:10	848-11	0000157	513	CT neg	GC POS	M
3/25/2011 18:37:10	848-12	0000158	751	CT neg	GC POS	M
3/25/2011 19:02:08	848-13	0000159	207	CT neg	GC POS	M
3/25/2011 19:02:09	848-14	0000160	811	CT neg	GC POS	M
3/25/2011 19:02:09	848-15	0000161	191	CT neg	GC POS	M
3/25/2011 19:02:09	947-01	0000162	558	CT neg	GC POS	M
3/25/2011 19:02:09	947-02	0000163	254	CT neg	GC POS	M
3/25/2011 19:12:08	947-03	0000164	1075	CT neg	GC POS	M
3/25/2011 19:12:09	947-04	0000165	436	CT neg	GC POS	M
3/25/2011 19:12:09	947-05	0000166	710	CT neg	GC POS	M
3/25/2011 19:12:09	947-06	0000167	679	CT neg	GC POS	M
3/25/2011 19:12:09	947-07	0000168	603	CT neg	GC POS	M
3/25/2011 19:17:08	947-08	0000169	759	CT neg	GC POS	M
3/25/2011 19:17:08	947-09	0000170	484	CT neg	GC POS	M
3/25/2011 19:17:09	947-10	0000171	1180	CT neg	GC POS	M
3/25/2011 19:17:09	947-11	0000172	1107	CT neg	GC POS	M

Printed: 3/28/2011 08:05:14

Page 6 of 7

Internal Training Use Only

CT/GC

Version: 1.33 (For In Vitro Diagnostic Use)

Time	Rack-Tube #	Sample ID	Total RLU	CT Result	GC Result	Flags
3/25/2011 19:17:09	947-12	0000173	1065	CT neg	GC POS	M