

RESULTS BY WORKLIST

Serial Number: 00239		System SW Version: 3.2.2.18		Date Range: 3/3/2012 12:53:24 PM to 3/3/2012 1:43:22 PM	
Assay: CT/GC					
Worklist IDs: 00239-20120303-01 Version: 1.36 (For In Vitro Diagnostic Use) ML: 606346 Exp: 5/15/2013					
Operators:					
Specimen Tested: 52		Invalid: 52			
Status Flag Legend: ci - Control Does Not Have Required Kinetic Profile. M - Sample ID of a new test request was manually entered via the keyboard when a barcode was not readable by the instrument. p - Assay Processing Error. x - Invalid Results Due To Control Failure.					
Reviewed By: _____		Date: _____		Verified By: _____ Date: _____	

CT/GC

Version: 1.36 (For In Vitro Diagnostic Use)

Worklist ID:

Time	Rack-Tube #	Sample ID	Total RLU	CT Result	GC Result	Flags
3/3/2012 12:53:24 PM	133-00	CT Positive Control - Lot 599663	672	CT POS	GC POS	ci
3/3/2012 12:53:24 PM	133-00	GC Positive Control - Lot 599392	705	CT neg	GC POS	
3/3/2012 12:53:24 PM	133-00	C12MC003379	35	Invalid	Invalid	x, p
3/3/2012 12:53:24 PM	133-00	C12MC003380	33	Invalid	Invalid	x, p
3/3/2012 12:53:24 PM	133-00	C12MC003381	29	Invalid	Invalid	x, p
3/3/2012 12:58:23 PM	133-00	C12MC003382	24	Invalid	Invalid	x, p
3/3/2012 12:58:23 PM	133-00	C12MC003383	38	Invalid	Invalid	x, p
3/3/2012 12:58:23 PM	133-00	C12MC003384	34	Invalid	Invalid	x, p
3/3/2012 12:58:23 PM	133-00	C12MC003385	32	Invalid	Invalid	x, p
3/3/2012 12:58:23 PM	133-00	C12MC003386	29	Invalid	Invalid	x, p
3/3/2012 1:03:23 PM	133-00	C12MC003387	43	Invalid	Invalid	x, p
3/3/2012 1:03:23 PM	133-00	C12MC003388	35	Invalid	Invalid	x, p
3/3/2012 1:03:23 PM	133-00	C12MC003389	39	Invalid	Invalid	x, p
3/3/2012 1:03:23 PM	133-00	C12MC003390	35	Invalid	Invalid	x, p
3/3/2012 1:03:23 PM	133-00	C12MC003391	34	Invalid	Invalid	x, p

CT/GC**Version: 1.36 (For In Vitro Diagnostic Use)**

Time	Rack-Tube #	Sample ID	Total RLU	CT Result	GC Result	Flags
3/3/2012 1:08:23 PM	146-00	C12MC003392	28	Invalid	Invalid	x, p
3/3/2012 1:08:23 PM	146-00	C12MC003393	32	Invalid	Invalid	x, p
3/3/2012 1:08:23 PM	146-00	C12MC003394	32	Invalid	Invalid	x, p
3/3/2012 1:08:23 PM	146-00	C12MC003395	26	Invalid	Invalid	x, p
3/3/2012 1:08:23 PM	146-00	C12MC003396	32	Invalid	Invalid	x, p
3/3/2012 1:13:23 PM	146-00	C12MC003397	27	Invalid	Invalid	x, p
3/3/2012 1:13:23 PM	146-00	C12MC003398	21	Invalid	Invalid	x, p
3/3/2012 1:13:23 PM	146-00	C12MC003399	23	Invalid	Invalid	x, p
3/3/2012 1:13:23 PM	146-00	C12MC003400	26	Invalid	Invalid	x, p
3/3/2012 1:13:23 PM	146-00	C12MC003401	19	Invalid	Invalid	x, p
3/3/2012 1:18:23 PM	146-00	C12MC003402	10	Invalid	Invalid	x, p
3/3/2012 1:18:23 PM	146-00	C12MC003403	13	Invalid	Invalid	x, p
3/3/2012 1:18:23 PM	146-00	C12MC003404	8	Invalid	Invalid	x, p
3/3/2012 1:18:23 PM	146-00	C12MC003405	1408	Invalid	Invalid	x, p
3/3/2012 1:18:23 PM	146-00	C12MC003406	9	Invalid	Invalid	x, p
3/3/2012 1:23:22 PM	030-00	C12MC003407	10	Invalid	Invalid	x, p
3/3/2012 1:23:22 PM	030-00	C12MC003408	8	Invalid	Invalid	x, p
3/3/2012 1:23:22 PM	030-00	C12MC003409	6	Invalid	Invalid	x, p
3/3/2012 1:23:22 PM	030-00	C12MC003410	7	Invalid	Invalid	x, p
3/3/2012 1:23:22 PM	030-00	C12MC003411	6	Invalid	Invalid	x, p
3/3/2012 1:28:22 PM	030-00	C12MC003412	17	Invalid	Invalid	x, p
3/3/2012 1:28:22 PM	030-00	C12MC003413	8	Invalid	Invalid	x, p
3/3/2012 1:28:22 PM	030-00	C12MC003414	7	Invalid	Invalid	x, p
3/3/2012 1:28:22 PM	030-00	C12MC003415	7	Invalid	Invalid	x, p
3/3/2012 1:28:22 PM	030-00	C12MC003416	1318	Invalid	Invalid	x, p
3/3/2012 1:33:22 PM	030-00	C12MC003417	21	Invalid	Invalid	x, p
3/3/2012 1:33:22 PM	030-00	C12MC003418	19	Invalid	Invalid	x, p
3/3/2012 1:33:22 PM	030-00	C12MC003419	9	Invalid	Invalid	x, p
3/3/2012 1:33:22 PM	030-00	C12MC003420	9	Invalid	Invalid	x, p
3/3/2012 1:33:22 PM	030-00	C12MC003421	5	Invalid	Invalid	x, p
3/3/2012 1:38:22 PM	267-00	C12MC003422	14	Invalid	Invalid	x, p
3/3/2012 1:38:22 PM	267-00	C12MC003423	5	Invalid	Invalid	x, p

CT/GC**Version: 1.36 (For In Vitro Diagnostic Use)**

Time	Rack-Tube #	Sample ID	Total RLU	CT Result	GC Result	Flags
3/3/2012 1:38:22 PM	267-00	C12MC003424	9	Invalid	Invalid	x, p
3/3/2012 1:38:22 PM	267-00	C12MC003425	10	Invalid	Invalid	x, p
3/3/2012 1:38:22 PM	267-00	C12MC003426	10	Invalid	Invalid	x, p
3/3/2012 1:43:22 PM	267-00	C12MC003427	46	Invalid	Invalid	x, p
3/3/2012 1:43:22 PM	267-00	C12MC003428	6	Invalid	Invalid	x, p
3/3/2012 1:43:22 PM	267-00	inhneg	8	Invalid	Invalid	M, x, p
3/3/2012 1:43:22 PM	267-00	inhpos	2555	Invalid	Invalid	M, x, p