

RESULTS BY WORKLIST

Serial Number: XXXXXXXX	System SW Version: 3.2.2.18	Date Range: 11/21/2011 1:23:13 PM to 11/21/2011 2:08:12 PM
Assay: HPV		
Worklist IDs: 00185-20111121-01 Version: 1.50 (For In Vitro Diagnostic Use) ML: 597942 Exp: 8/15/2012 Analyte Cutoff: 0 IC Cutoff: 0		
Operators: XXXXXXXX		
Specimen Tested: 40	Invalid: 40	
Status Flag Legend: f - Calibration failed. p - Assay processing error.		
Reviewed By: _____	Date: _____	Verified By: _____ Date: _____

HPV

Version: 1.50 (For In Vitro Diagnostic Use)

Worklist ID: 00185-20111121-01

Time	Rack-Tube #	Sample ID	IC RLU	IC Result	RLU	S/CO	Result	Flags
11/21/2011 1:23:13 PM	574-01	Negative Calibrator - Lot 588262	269203	Invalid	0	0.00	Invalid	p
11/21/2011 1:23:13 PM	574-01	Negative Calibrator - Lot 588262	254375	Invalid	0	0.00	Invalid	p
11/21/2011 1:23:13 PM	574-01	Negative Calibrator - Lot 588262	251421	Invalid	0	0.00	Invalid	p
11/21/2011 1:23:13 PM	574-02	Positive Calibrator - Lot 588263	248574	Invalid	1296525	0.00	Invalid	p
11/21/2011 1:23:13 PM	574-02	Positive Calibrator - Lot 588263	247074	Invalid	1288704	0.00	Invalid	p
11/21/2011 1:28:13 PM	574-02	Positive Calibrator - Lot 588263	279351	Invalid	1321475	0.00	Invalid	p
11/21/2011 1:28:13 PM	574-03	HPV-PANEL-A-591503	264351	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:28:13 PM	574-04	HPV-PANEL-A-591503	240310	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:28:13 PM	574-05	HPV-PANEL-A-591503	249805	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:28:13 PM	574-06	HPV-PANEL-A-591503	252168	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:33:13 PM	574-07	HPV-PANEL-A-591503	271026	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:33:13 PM	574-08	HPV-PANEL-A-591503	255852	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:33:13 PM	574-09	HPV-PANEL-A-591503	255594	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:33:13 PM	574-10	HPV-PANEL-A-591503	262209	Invalid	0	0.00	Invalid	f, p

HPV**Version: 1.50 (For In Vitro Diagnostic Use)**

Time	Rack-Tube #	Sample ID	IC RLU	IC Result	RLU	S/CO	Result	Flags
11/21/2011 1:33:13 PM	574-11	HPV-PANEL-A-591503	247759	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:38:13 PM	574-12	HPV-PANEL-A-591503	262731	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:38:13 PM	574-13	HPV-PANEL-A-591503	247904	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:38:13 PM	574-14	HPV-PANEL-A-591503	245938	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:38:13 PM	574-15	HPV-PANEL-A-591503	249686	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:38:13 PM	680-01	HPV-PANEL-A-591503	253718	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:43:12 PM	680-02	HPV-PANEL-A-591503	272139	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:43:12 PM	680-03	HPV-PANEL-A-591503	256261	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:43:12 PM	680-04	HPV-PANEL-A-591503	242139	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:43:12 PM	680-05	HPV-PANEL-A-591503	259426	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:43:13 PM	680-06	HPV-PANEL-A-591503	248911	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:48:12 PM	680-07	HPV-PANEL-A-591503	283507	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:48:12 PM	680-08	HPV-PANEL-A-591503	261720	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:48:12 PM	680-09	HPV-PANEL-A-591503	243055	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:48:12 PM	680-10	HPV-PANEL-A-591503	250646	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:48:12 PM	680-11	HPV-PANEL-A-591503	252430	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:53:12 PM	680-12	HPV-PANEL-A-591503	263492	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:53:12 PM	680-13	HPV-PANEL-A-591503	255321	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:53:12 PM	680-14	HPV-PANEL-A-591503	243181	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:53:12 PM	680-15	HPV-PANEL-A-591503	248244	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:53:12 PM	457-01	HPV-PANEL-A-591503	241806	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:58:12 PM	457-02	HPV-PANEL-A-591503	276765	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:58:12 PM	457-03	HPV-PANEL-A-591503	259846	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:58:12 PM	457-04	HPV-PANEL-A-591503	253564	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:58:12 PM	457-05	HPV-PANEL-A-591503	250841	Invalid	0	0.00	Invalid	f, p
11/21/2011 1:58:12 PM	457-06	HPV-PANEL-A-591503	255297	Invalid	0	0.00	Invalid	f, p
11/21/2011 2:03:12 PM	457-07	HPV-PANEL-A-591503	246725	Invalid	0	0.00	Invalid	f, p
11/21/2011 2:03:12 PM	457-08	HPV-PANEL-A-591503	244869	Invalid	0	0.00	Invalid	f, p
11/21/2011 2:03:12 PM	457-09	HPV-PANEL-A-591503	252313	Invalid	0	0.00	Invalid	f, p
11/21/2011 2:03:12 PM	457-10	HPV-PANEL-A-591503	246848	Invalid	0	0.00	Invalid	f, p
11/21/2011 2:03:12 PM	457-11	HPV-PANEL-A-591503	243852	Invalid	0	0.00	Invalid	f, p
11/21/2011 2:08:12 PM	457-12	HPV-PANEL-A-591503	252715	Invalid	0	0.00	Invalid	f, p

HPV

Version: 1.50 (For In Vitro Diagnostic Use)

For Internal Use Only