



QIASTAT-DX (QIAGEN)
Preventive Maintenance
CHECK LIST



Customer Name: _____

Operational Module SN: _____

Analytical Modules SN: _____

No.	Operation	Action
1.	Cleaning & Disinfection	
	Clean the indicated surfaces (Using 70 % ethanol /suitable cleaning agent) and damp cloth	Instrument surfaces
		OM. Touchscreen
		Barcode window
2.	Software/Assay	Check the software/firmware version
		Check the assays version
		Update if required
3.	Air Filter Exchange for AMs	Cartridge exchange
4.	Operational check and verification by Infield Service Tool for each AM	Visual inspection
		Actuators Status
		Pneumatics
		qPCR internal standard
		Thermal unit & Liquid sensors test
5.	Testing	Test the Modules by running ACC if applicable
6.	Final testing	Request the user to perform run using patient samples

Date: _____

Engineer Name: _____

Signature: _____