

GRIFOLS Diagnostic Grifols, S.A.	Reference:	REGD-0000790	Version:	4.0, CURRENT
	Status:	Effective	Effective Date:	26-May-2011
Title:	SM DG57-DG517. Serial number configuration.			Page 2 of 4

DG-57 / DG-517 Service Manual

SERIAL NUMBER CONFIGURATION

Procedure description
Every instrument is identified by a serial number 057-NNNNNNNN in case of DG-57 and 517-NNNNNNNN in case of DG-517. Serial number is located in the rate plate (the label placed in the back of the instrument). NNNNNNN is a sequential number. In order to know what instrument is connected to the PC, the first four 'N' digits of the serial number must be stored in the Carousels card EEPROM. This number is then used by the software to identify batch and log files generated by this instrument. This number is saved during manufacturing process, and under normal circumstances it must never be saved again. However, if a EEPROM failure occurs and it must be replaced, or the Carousels card is replaced, the serial number must be saved again. This procedure explains how to save the serial number of an instrument. <i>Note: Document REGD-0000790_v4 replaces document REGD-0000790_v3.</i>

GRIFOLS Diagnostic Grifols, S.A.	Reference:	REGD-0000790	Version:	4.0, CURRENT
	Status:	Effective	Effective Date:	26-May-2011
Title:	SM DG57-DG517. Serial number configuration.			Page 3 of 4

EQUIPMENT

1. An instrument connected with a PC with the software installed

SET UP

1. Turn on the instrument.
2. Execute the ST software (Technical Service software).
3. Enter the Technical Service password (d+g_diana).

PROCEDURE

Process	Expected Results
1. Locate instrument's serial number in the rate plate at the back of the instrument (057-NNNNNNN or 517-NNNNNNN).	
2. Access 'Settings' tab in ST software.	---
3. Locate the "Ser N°" frame in the lower right part of the 'Settings' tab.	---
4. Type instrument's serial number first four 'N' digits in the "Ser N°" textbox.	The number can be changed.
5. Exit the software pressing the "Exit" button.	---

GRIFOLS Diagnostic Grifols, S.A.	Reference:	REGD-0000790	Version:	4.0, CURRENT
	Status:	Effective	Effective Date:	26-May-2011
Title:	SM DG57-DG517. Serial number configuration.			Page 4 of 4

DOCUMENT CHANGE HISTORY

Document change history			
Version	Date	Change	Rationale
2.0	2005/05/06	<i>Document REGD-0000790_v2 replaces document 57018V0M000.</i>	Document reference changed
3.0	2011/05/02	Added `previous version` `data`, `change` and `rationale` inside document change history section.	Update change history for previous document history modification.
3.0	2011/05/02	Updated the Document Change History: Date format changed to YYYY/MM/DD Dates changed to document version effective date.	Adoption of ISO 8601 date format to avoid ambiguity.
3.0	2011/05/02	ID-GELSTATION and PROVUE erased from the title.	Service manual intended for Grifols market.
3.0	2011/05/02	<i>Document REGD-0000790_v3 replaces document REGD-0000790_v2.</i>	Document version changed.
3.0	2011/05/02	Updated serial number format inside `procedure description` section.	Serial number format changed.
3.0	2011/05/02	Update `procedure` section due to the new serial number format.	Serial number format changed.
4.0	See document header	Document title.	Commercial reference substituted by project reference.
4.0	See document header	<i>Document REGD-0000790_v4 replaces document REGD-0000790_v3.</i>	Document version changed.
4.0	See document header	Procedure description section and procedure section updated.	Document available for DG-57 and DG-517.
4.0	See document header	Header title.	Increase information to easy document search.