

GRIFOLS Diagnostic Grifols, S.A.	Reference:	REGD-0000802	Version:	3.0, CURRENT
	Status:	Effective	Effective Date:	16-Jun-2011
Title:	SM DG57-DG517. PCA centrifuge-2 control current limits adjustment.			Page 2 of 4

DG-57 / DG-517 Service Manual

PCA CENTRIFUGE-2 CONTROL CURRENT LIMITS ADJUSTMENT

Procedure description
This procedure describes how adjust the electronics values of Centrifuge-2 PCA. <i>Note: Document REGD-0000802_v3 replaces document REGD-0000802_v2.</i>

GRIFOLS Diagnostic Grifols, S.A.	Reference:	REGD-0000802	Version:	3.0, CURRENT
	Status:	Effective	Effective Date:	16-Jun-2011
Title:	SM DG57-DG517. PCA centrifuge-2 control current limits adjustment.			Page 3 of 4

EQUIPMENT

- Multimeter.



NOTE:

Ground in all measurements is in TP20. (see Picture 1)

SET UP

- Equipment turned on.

PROCEDURE

Process	Expected Results
Adjusting X-axis current limit - Measure the voltage in TP24 with the multimeter (see Picture 1). - If required, adjust P7 until the measured value reaches expected result (see Picture 1).	990 ± 10 mV
Picture 1	990 ± 10 mV

GRIFOLS Diagnostic Grifols, S.A.	Reference:	REGD-0000802	Version:	3.0, CURRENT
	Status:	Effective	Effective Date:	16-Jun-2011
Title:	SM DG57-DG517. PCA centrifuge-2 control current limits adjustment.			Page 4 of 4

DOCUMENT CHANGE HISTORY

Document change history			
Version	Date	Change	Rationale
2.0	2005/05/09	Document REGD-0000802_v2 replaces document 57030V0M000.	Document reference changed.
3.0	See document header.	Added previous version 'data', 'change' and 'rationale' inside document change history section.	Update change history for previous document history modification.
3.0	See document header.	Document REGD-0000802_v3 replaces document REGD-0000802_v2.	Document version changed.
3.0	See document header.	Document title.	Commercial reference substituted by project reference.
3.0	See document header.	Document header.	Increase information to easy document search.