

GRIFOLS Diagnostic Grifols, S.A.	Reference:	REGD-0000838	Version:	4.0, CURRENT
	Status:	Effective	Effective Date:	23-May-2011
Title:	SM DG57-DG517. Files and folders structure.			Page 2 of 4

DG-57 / DG-517 Service Manual

FILES AND FOLDERS STRUCTURE

Procedure description
This document describes the files and folders structure of the software. The default installation of the Diana System Software creates a root directory with the program and definition files, and nine subdirectories.

GRIFOLS Diagnostic Grifols, S.A.	Reference:	REGD-0000838	Version:	4.0, CURRENT
	Status:	Effective	Effective Date:	23-May-2011
Title:	SM DG57-DG517. Files and folders structure.			Page 3 of 4

EQUIPMENT

1. ---

SET UP

1. ---

PROCEDURE

Process	Expected Results
Contents of the root Directory 1. Applications: • WADiana.exe: DG-57 / DG-517 software • Lector.exe: Diana R software • Listas.exe: Diana O software • Config.exe: configuration program • DianaServer.exe: connection to Host program • Etiquetes.exe: labels generator • BCGen.exe: labels with format generator • WADDiagnostic.exe: WADDiagnostic tool • VisioServer.exe: image processing internal routines Additionally HostDemo.exe can be found inside software CD: • HostDemo.exe: Host simulator. Hostdemo will be available in the CD path 'Programs\Server'. To install it just copy HostDemo.exe file to Diana System installation folder. Definition files: • Debug.def: image procession parameters and options • Idioma*.def: multilanguage files for all applications Other files, *.dll, *.bpl, ...	---
Folders 2. CQC folder: • Quality Control files	---
3. Log folder: • WADiana1.log: log of WADiana C Software • WADiana1.old: backup of WADiana1.log file when it exceeds 1Mb • DianaO.log: log of Diana O Software • DianaR.log: log of Diana R Software • VisioServer.log: log of VisioServer Software • Tolcent.txt: position control of the centrifuge tolerance statistics • Debug.txt: contains the last 100 instructions sent through the COM port	---
4. Defs folder: • *.def: Definition files for cards, reagents, interpretation of results,...	---
5. Grus folder: • *.gru: Test definition files	---

GRIFOLS Diagnostic Grifols, S.A.	Reference:	REGD-0000838	Version:	4.0, CURRENT
	Status:	Effective	Effective Date:	23-May-2011
Title:	SM DG57-DG517. Files and folders structure.			Page 4 of 4

Process		Expected Results
6.	Host folder: Host.db: host transmission database	---
7.	Plantillas folder: *.pln: templates definition files	---
8.	Tandas: *.tan: files with the information of a batch *.bmp: wells pictures of all cards. Their name is the combination of the batch number, card number and well number	---
9.	Usades: Tarjetas.db: database with partially used cards so they can be reused in other batches.	---
10.	Tars: *.tar: test definitions for manual batches created with Diana Reader.	---

DOCUMENT CHANGE HISTORY

Document change history			
Version	Date	Change	Rationale
2.0	2005/05/09	Document REGD-0000838_v2 replaces document 57069V0M000.	Uploaded document to GxPharma.
3.0	2011/04/28	Eliminated the reference to ID-GelStation and ProVue distributions.	Service manual intended for Grifols market'.
3.0	2011/04/28	Updated the Document Change History: - Date format changed to YYYY/MM/DD. - Dates changed to document version effective date.	Adoption of ISO 8601 date format to avoid ambiguity.
3.0	2011/04/28	HostDemo.exe now is not installed from CD room to root directory. It will remain inside CD room.	Software erased from standard installation due to that it is a development and verification tool.
4.0	See document header	Document REGD-0000838_V4 replaces REGD-0000838_V3.	Document version changed.
4.0	See document header	-Document title. -Point 1 'applications' description modified.	Commercial reference substituted by project reference.
4.0	See document header		
4.0	See document header	Header title.	Increase information to easy document search.