

GRIFOLS Diagnostic Grifols, S.A.	Reference:	REGD-0000847	Version:	3.0, CURRENT
	Status:	Effective	Effective Date:	26-May-2011
Title:	SM DG57-DG517. Log file description.			Page 2 of 9

DG-57 / DG-517 Service Manual

LOG FILE DESCRIPTION

Procedure description
This procedure describes the log file and how to read it. <i>Note: Document REGD-0000847_v3 replaces document REGD-0000847_v2.</i>

GRIFOLS Diagnostic Grifols, S.A.	Reference:	REGD-0000847	Version:	3.0, CURRENT
	Status:	Effective	Effective Date:	26-May-2011
Title:	SM DG57-DG517. Log file description.			Page 3 of 9

EQUIPMENT

- | |
|-------------------------------------|
| 1. A PC with the software installed |
|-------------------------------------|

SET UP

- | |
|-------------------|
| 1. Turn on the PC |
|-------------------|

PROCEDURE 1 – Locate and open the Log file

Process	Expected Results
1. Locate the installation directory.	---
2. Open \Log directory under the installation directory.	---
3. Open the file wadiana1.log with Notepad.	---

PROCEDURE 2 – Interpreting the Log file

Log entry description	Comment
1. Log file is electronically signed. First entry in log file can be ## LOG SIGNATURE OK!! ## or ## LOG SIGNATURE FAILED!! ##	## LOG SIGNATURE FAILED!! ## appears when the log file doesn't exist previously, the log file wasn't signed or the log file was modified. ---
2. The following entries are [I] Free Space : LOCAL_DRIVE: 24% (3561728 KB) ## START LOG ##	---
3. Next, the software starts to log the initialization of the instrument *** Begin WADiana ver. X.XX *** Initializing IO Initializing Cavro plunger Initializing XYZ Initializing Centrifuge Initializing Carrousel Reading IO configurations *** TAULA Calibration values *** Empty Waste : 417 Full Waste : 2275 Empty A : 396 ...	In this part of log file we can identify information like the software version, instrument serial number and calibration values.
4. After initializing, incubator temperatures are set to default values and the user can start to place orders. Adjusting Temperatures: (Temp1: 25 - Temp2: 25) Door Closed User: START Temperatures adjusted (Temp1: 24.2 - Temp2: 24.4) User:Manual tube ID New tube in position 1, ID: , Size: 0	Most user actions starts with "User:"

GRIFOLS Diagnostic Grifols, S.A.	Reference:	REGD-0000847	Version:	3.0, CURRENT
	Status:	Effective	Effective Date:	26-May-2011
Title:	SM DG57-DG517. Log file description.			Page 4 of 9

Log entry description	Comment
5. When a user starts a worksheet samples, reagents and card identification are performed. <pre>User: START TestWiz ***** Begin worksheet ***** WorkSheet launched by user (null) Begin samples identification Identified sample 1975270463 in position 1, 15 width End samples identification Begin reagents identification Begin Cards identification Card 0135020211104206 identified on position 1 End Cards identification Vial detected in position 9 End reagents identification New batch 372. GRU: Wad02A.gru New patient :1975270463 . Serum in position 1 with diameter 15 , Red cells in position 1 with diameter 15 . New reagent 0310(M) in position 9 Begin primming End primming</pre>	---
6. Each pipeting block is described. <pre>Begin dilution for batch 372 Aspiring Reactivo DIL 1 en la posición 9 Aspiring Hematíes de la muestra 1 Dispensing to dilutions station End dilution Begin Pipeting for batch 372 Aspiring Dilución de la muestra 1 Dispensing Dilución de la muestra 1 Well A.GD, 1 in 0135020211104206 Well B.GD, 2 in 0135020211104206 Well D.RH, 3 in 0135020211104206 Well ctl.GD, 4 in 0135020211104206 End Pipeteing</pre>	---
7. Incubation and centrifugation are logged. <pre>Begin waiting for incubation of 1 cards for batch 372, subbatch 1 Begin Dilution Station washing Ready for centrifuging 1 cards of batch 372, subbatch 1 Temp1: 24.9 - Temp2: 24.8 End incubation of 1 cards for batch 372, subbatch 1 End Dilution Station washing Loading 0135020211104206 card into the centrifuge Loading balancing card into the centrifuge Begin centrifuging 1 cards of batch 372, subbatch 1 Start Centrifuging at 1840 RPMs during 540 seconds (nominal values) Start Centrifuging at 1843.1 RPMs during 553.0 seconds (real values) Centrifuged at 1839.3 RPMs End centrifuging 1 cards of batch 372, subbatch 1</pre>	Temperatures and time and speed of centrifugation are shown. "Nominal values" means those than appear in the test file (*.GRU) and "real values" are those really programmed by de software due to the electronics resolution. Finally an estimated value by the software of the real centrifugation speed is shown
8. After centrifuging, the cards are read and the results are shown. <pre>Begin Reading Readed Card : 0135020211104206 * Results : 1+ 1+ 1+ 2+ * Volumes : 50 51 53 47 50 0 0 0 Attention : Doubtful card in 1 End Reading User: START , Cards Incidences ... User: Open door Door Closed Doubtful resolution 0 : Some wells are unused User: AGREE TestWiz</pre>	Estimation of the dispensed volume in each well is displayed. Unused wells have no result and 0 volume.
9. When the software closes, ends the log file with the last entry and the digital signature. <pre>## END LOG ## 3C266CB5954F5016821C182CA428DF1D8C94BA73184CA63B1D60B5CB293...</pre>	---

GRIFOLS Diagnostic Grifols, S.A.	Reference:	REGD-0000847	Version:	3.0, CURRENT
	Status:	Effective	Effective Date:	26-May-2011
Title:	SM DG57-DG517. Log file description.			Page 5 of 9

Annex A – Alphabetically sorted list of all log entries

Log Entry	Description	Type
**** BRAZO Calibration values ****	Calibrations values stored by BRAZO module will follow.	Execution trail
*** Begin #Software Name# ver. #Software version# ***	The Wadiana software begins	Execution trail
*** CARRU Calibration values ***	Calibrations values stored by CARRUSEL module will follow.	Execution trail
*** End #Software name#***	The Wadiana software ends in a controlled manner	Execution trail
*** TAULA Calibration values ***	Calibrations values stored by TAULA module will follow.	Execution trail
*** VISIO Calibration values ***	Calibrations values stored by VISIO module will follow.	Execution trail
***** Begin worksheet *****	A new worksheet begins execution.	Execution trail
Abort cards cancelling	A cancelling cards process aborts.	Errors trail
Abort centrifuging	A centrifuging process is aborted	Errors trail
Abort centrifuging %d cards of batch %d, subbatch %d	Subbatch centrifuging process aborted.	Errors trail
Abort Execution	The execution of a subbatch was aborted.	Errors trail
Abort priming	Priming procedure has been aborted	Execution trail
Abort Reading	A reading process aborts.	Execution trail
Abort Test	A "stop" button in an error red window was pressed and the test will be aborted.	Errors trail
Aborting centrifugation due open door and the arm don't block centrifuge	The door was detected to be open during a centrifuge process and the arm isn't blocking the centrifuge entrance.	Errors trail
Aborting incubation of %d cards for batch %d, subbatch %d	An incubation process is aborted.	Execution trail
Adjusting Temperatures: (Temp1: %d - Temp2 %d)	Temperatures aren't in the right range yet.	Execution trail
Aspiring %s	A sample, reagent or dilution will be aspired.	Execution trail
Attention : Doubtful card in %d	A doubtful card was placed in the service rack.	Execution trail
B1 24 : %ld	Calibration for incubator bloc 1 at 24 °C	Execution trail
B1 37 : %ld	Calibration for incubator bloc 1 at 37 °C	Execution trail
B2 24 : %ld	Calibration for incubator bloc 2 at 24 °C	Execution trail
B2 37 : %ld	Calibration for incubator bloc 2 at 37 °C	Execution trail
Base : %lu	Calibration for the centrifuge positioning.	Execution trail
Begin Cards identification	A cards identification process begins.	Execution trail
Begin centrifuge loop	Technical service "Load Centrifuge" procedure begins.	Execution trail
Begin centrifuging %d cards of batch %d, subbatch %d	A centrifuge process will begin for a subbatch.	Execution trail
Begin dilution for batch %d	A dilution block begins.	Execution trail
Begin Dilution Station washing	A dilution station washing procedure starts.	Execution trail
Begin Pipeting for batch %d	A pipeting block begins.	Execution trail
Begin priming	Priming procedure had started.	Execution trail
Begin Reading	A reading process begins.	Execution trail
Begin reagents identification	A reagents identification process begins.	Execution trail
Begin samples identification	A sample identification process begins.	Execution trail
Begin waiting for incubation of %d cards for batch %d, subbatch %d	Starts a wait for incubation process.	Execution trail
Bright cam%d : %ld	Brightness calibration of a camera.	Execution trail
Calibration for holder %d: %d	Absolute calibration for a holder during centrifuge calibration	Execution trail
Calibration offset for holder %d: %d	Offset calibration for a holder during centrifuge calibration	Execution trail
Can not release the card	A card can't be released from the gripper.	Errors trail
Cancel by clot	A sample is cancelled by clot.	Errors trail
Cancel Diluent!	Diluent cancelled during a dilution block.	Errors trail
Cancel Sample!	Sample cancelled during a dilution block.	Errors trail
Cancelling %s. Error ???	A sample is cancelled for an unexpected error.	Errors trail
Cancelling %s. Level detection fail	A sample or reagent was cancelled due to level detection fail.	Errors trail
Cancelling %s. Tube not found	A tube was removed from the carrousel before the execution is finished.	Errors trail
Cancelling well	A well was cancelled.	Errors trail
Can't read the card	The program can't read a card during the identification.	Errors trail
Can't throw the card to the trash	The gripper can't release the card to throw it to the trash	Errors trail
Card %s identified on position %d	A card was identified. The position of the card in the incubator is shown.	Execution trail
Card in wrong heating block	The program has rejected a card during the identification because it is the wrong heating block.	Errors trail
Cards cancelling	A centrifuge process was cancelled and all cards in centrifuge are cancelled.	Execution trail
CarrErrCod %c	An error code was returned by the "Carrousel" card. See the technical manual for a error description.	Errors trail
CentErrCod %c	An error code was returned by the "Centrifuge" card. See the technical manual for a error description.	Errors trail
Centrifugation velocity out of range: %6.1lf RPMS (%6.1lf - %6.1lf)	A centrifugation procedure has ended unsuccessfully. Valid RPMS range is logged.	Errors trail
Centrifuge calibration procedure	Centrifuge calibration procedure has started	Execution trail
Centrifuge loaded	Technical service "Load Centrifuge" procedure has loaded the centrifuge.	Execution trail
Centrifuge loop # %d	Technical service "Load Centrifuge" procedure loop counter.	Execution trail
Centrifuge retries positioning	The centrifuge retries positioning in order to get more accuracy.	Execution trail

GRIFOLS Diagnostic Grifols, S.A.	Reference:	REGD-0000847	Version:	3.0, CURRENT
	Status:	Effective	Effective Date:	26-May-2011
Title:	SM DG57-DG517. Log file description.			Page 6 of 9

Log Entry	Description	Type
Centrifuge unloaded	Technical service "Load Centrifuge" procedure has unloaded the centrifuge.	Execution trail
Centrifuged at %6.1lf RPMs	A centrifugation procedure has ended successfully	Execution trail
Clot detected during a dilution process	A clot was detected during a dilution process	Errors trail
Clot in Serum	A clot is detected in a serum tube.	Errors trail
Code doesn't match	The program has rejected a card during the identification because bar code doesn't match	Errors trail
ComErrCod %d	An error code was returned by the "Communications" module. See the technical manual for a error description.	Errors trail
Contrast cam%d : %ld	Contrast calibration of a camera.	Execution trail
Dilutions station clogged! Aborting!	The dilutions station was clogged due a clot.	Errors trail
Dispensing %s	A sample, reagent or dilution will be dispensed.	Execution trail
Dispensing to dilutions station	The dilution is dispensed to the dilutions station	Execution trail
Door Closed	The door is closed.	Execution trail
Doubtful resolution %d : %s	The user acknowledged a card placed in the doubtful rack. The card position and the cause are given.	User action
Empty A : %ld	Calibration for empty A bottle.	Execution trail
Empty B : %ld	Calibration for empty B bottle.	Execution trail
Empty Waste : %ld	Calibration for empty Waste bottle.	Execution trail
End cards cancelling	A cancelling cards process ends.	Execution trail
End Cards identification	A cards identification process ends.	Execution trail
End centrifuge loop	Technical service "Load Centrifuge" procedure ends.	Execution trail
End centrifuging %d cards of batch %d, subbatch %d	A centrifuge process was ended for a subbatch.	Execution trail
End dilution	A dilution block was ended.	Execution trail
End Dilution Station washing	A dilution station washing procedure ends	Execution trail
End incubation of %d cards for batch %d, subbatch %d	An incubation process ends.	Execution trail
End Pipeteing	A pipeting block was ended.	Execution trail
End priming	Priming procedure had ended.	Execution trail
End Reading	A reading process ends.	Execution trail
End reagents identification	A reagents identification process ends.	Execution trail
End samples identification	A sample identification process ends.	Execution trail
End Tubes Identification	The samples identification process started by user has ended.	User action
Error : Can't handle the card	A card was detected when the gripper try to get it, but it can't handle it.	Errors trail
Error : Can't left the card	A handled card by the gripper can't be left.	Errors trail
Error : Card not found	A card wasn't detected when the gripper try to get it.	Errors trail
Error : Lost card	The gripper doesn't detect any more a handled card.	Errors trail
Error Diluent : %s	Error detected when aspiring Diluent. See the technical manual for a error description.	Errors trail
Error during centrifugation : %s	An error was occurred during centrifugation. See the technical manual for a error description.	Errors trail
Error Reader. Can't read the card. %s	A card can't be read. See the technical manual for a error description.	Errors trail
Error Red cells : %s	Error detected when aspiring red cells. See the technical manual for a error description.	Errors trail
Expired Card	The program has rejected a card during the identification because it is expired.	Errors trail
Full A : %ld	Calibration for full A bottle.	Execution trail
Full B : %ld	Calibration for full B bottle.	Execution trail
Full Waste : %ld	Calibration for full Waste bottle.	Execution trail
General exception %s	Error not expected	Errors trail
Holder #holder number# offset: %ld	Calibration relative to each centrifuge holder.	Execution trail
Identified sample %s in position %d, %d width	A sample was identified in the shown position of the samples carrousel. The bar code and the tube width are also shown	Execution trail
Identified vial %s in position %d	A reagent was identified in the shown position of the reagents carrousel. The bar code is shown.	Execution trail
Initializing Carrousel	Carrousel card will initialise.	Hardware trail
Initializing Cavro plunger	Cavro plunger will initialise.	Hardware trail
Initializing Centrifuge	Centrifuge card will initialise.	Hardware trail
Initializing IO	IO card will initialise.	Hardware trail
Initializing XYZ	XYZ card will initialise.	Hardware trail
Instrument initialized. Ser. N° %ld	The instrument was successfully initialised. Shows serial number of the instrument.	Hardware trail
IOErrCod %c	An error code was returned by the "IO" card. See the technical manual for a error description.	Errors trail
LecErrCod : %d	An error code was returned by the "Reader" module. See the technical manual for a error description.	Errors trail
Loading %s card into the centrifuge	A card is loaded into the centrifuge.	Execution trail
Loading balancing card into the centrifuge	The balancing card is loaded into the centrifuge.	Execution trail
Magnetic shaker: %s	There are magnetic shaker?	Execution trail
MaqErrCod %d	An error code was returned by the "Maquina" module. See the technical manual for a error description.	Errors trail
Missing %d cards of type %s	There are not enough cards in order to start the execution.	Errors trail
New batch %d. GRU: %s	A new batch begins. The batch number and the GRU file name are shown.	Execution trail

GRIFOLS Diagnostic Grifols, S.A.	Reference:	REGD-0000847	Version:	3.0, CURRENT
	Status:	Effective	Effective Date:	26-May-2011
Title:	SM DG57-DG517. Log file description.			Page 7 of 9

Log Entry	Description	Type
New donor :%s. Serum in position %d %s with diameter %d %s, Red cells in position %d %s with diameter %d %s.	A new donor is included in the new batch. Tube positions and widths are shown.	Execution trail
New patient :%s. Serum in position %d %s with diameter %2.0f %s, Red cells in position %d %s with diameter %2.0f %s.	A new patient is included in the new batch. Tube positions and widths are shown.	Execution trail
New reagent %s %s in position %d	A reagent is used in the new batch. The barcode and the position is shown.	Execution trail
New tube in position #tube position#, ID: #Barcode# #Manually ID?#, Size: #tube diameter#	A new tube had been identified.	User action
No existent tube in position #Tube position#	The tube in that position is no longer there and user will be asked to delete it.	Execution trail
No level in %s	There was no enough level of liquid for pipeting.	Errors trail
Not enough volume for pipeting all wells	There are not enough liquid to pipet all wells. Less wells will be pipeted in order to optimize liquid.	Execution trail
Over Incubation!	The incubator went out of temperature range.	Errors trail
Over Temperature!!!!!!	A over temperature error has detected.	Errors trail
Query Error : #Error code#	A error red window was pop-up showing an error. See Technical Manual for a description of all errors.	Execution trail
Read Card : %s * Results :	A card is read and shows the results and the dispensed volumes detected by the reader.	Execution trail
Reading Centrifuge configurations	The EEPROM of the Centrifuge card will be read.	Hardware trail
Reading IO configurations	The EEPROM of the IO card will be read.	Hardware trail
Reading Vision configurations	Visio parameters will be read from the EEPROM.	Hardware trail
Reading XYZ configurations	The EEPROM of the XYZ card will be read.	Hardware trail
Ready for centrifuging %d cards of batch %d, subbatch %d	A subbatch pipeting and incubation was finished and all its cards are ready to centrifuge.	Execution trail
Reagents offset : %ld	Offset for the reagents carousel	Execution trail
Sample Camera distance : %d	Distance of the sample's camera.	Execution trail
Sample cancelled by clot	A sample was cancelled due a clot.	Execution trail
Samples CB -- #Readings: %ld corrects: %ld NoCodeReaded: %ld Errors: %ld	Technical service "Samples reading" procedure results.	Execution trail
Serial num. : %d	Serial number of the instrument.	Execution trail
Some wells aren't empty	The program has rejected a card during the identification because it some wells aren't empty	Errors trail
Start Centrifuging at %d RPMs during %d seconds (nominal values)	A centrifugation procedure starts. Nominal values are logged.	Execution trail
Start Centrifuging at %d RPMs during %d seconds (real values)	A centrifugation procedure starts. Real values are logged.	Execution trail
Start incubation of %d cards for batch %d, subbatch %d	A incubation starts.	Execution trail
Start timer in batch %d, subbatch %d	A timer starts to count in a batch	Execution trail
STAT time out. Beep sounds.	An STAT pause was very long.	Execution trail
Temp1: %3.1lf - Temp2: %3.1lf	Measured temperatures of the incubation blocks for a subbatch.	Execution trail
Temperature out of range (Temp1: %3.1lf - Temp2: %3.1lf)	Temperatures goes out of range during incubation.	Errors trail
Temperature time out (Temp1: %3.1lf - Temp2: %3.1lf)	Temperatures didn't reach the right values in the expected time.	Errors trail
Temperatures adjusted (Temp1: %3.1lf - Temp2: %3.1lf)	Temperatures had reached the right values.	Execution trail
There is no sample in position %d	A sample was missing during the identification in the shown position of the samples carousel	Execution trail
There is no tube in position %d	A tube was missing in the shown position of the samples carousel	Execution trail
There is no vial in position %d	A vial was missing in the shown position of the reagents carousel	Execution trail
Timer ends in batch %d, subbatch %d	A timer ends and the execution continues.	Execution trail
Too much short centrifugation to measure velocity	The centrifugation is too much short to measure a precise velocity. Typically shown in technical service routines.	Execution trail
Tube confirmed in position %d	A tube was detected in the shown position of the samples carousel	Execution trail
Tube deleted in position #Tube position#	The user had deleted a non existent tube.	User action
Tube in position #tube position#, ID: #Barcode# #Manually ID?#, Size: #tube diameter#	A tube was confirmed or modified by the user.	User action
Unhandled %s	Error not expected.	Errors trail
User #Username# Logout	The user #Username# has chosen the "Logout" option of the pop-up menu.	User action
User configuration Log : #Username#	The user #Username# has validated his password and had permission enter in the configuration window.	User action
User Log : %s	The user enters a correct password and log-in.	User action
User Log fail	The user has entered a wrong password or hadn't sufficient privileges.	User action
User QC Log : #Username#	The user #Username# has validated his password and had permission enter in the QC window.	User action
User: #tube diameter# tube selected	The user had sleeted to fix the tube diameter to #tube diameter#	User action
User: Abort TestWiz Correction	The user pushed the "cross" button in the selection test window during a identification incidence report.	User action
User: Abort TestWiz Incidences	The user pushed the "cross" button in the selection test window during a final incidences report.	User action
User: Abort TestWiz New test	The user pushed the "cross" button in the selection test window.	User action

GRIFOLS Diagnostic Grifols, S.A.	Reference:	REGD-0000847	Version:	3.0, CURRENT
	Status:	Effective	Effective Date:	26-May-2011
Title:	SM DG57-DG517. Log file description.			Page 8 of 9

Log Entry	Description	Type
User: Abort Urgency	The user pushed the "cross" button in the selection test window during a urgency.	User action
User: AGREE TestWiz	The user agree a incidence reported about the just finished execution.	User action
User: Automatic tube detection selected	The user had selected automatic tub diameter detection	User action
User: Cancel Emergency open and continue	User has cancelled an emergency open while waiting for open the door.	User action
User: Configuration	The user has chosen the "Configuration" option of the pop-up menu.	User action
User: Continue	The user pushed the "arrow" button in the stop execution window.	User action
User: Emergency Open	The user pushed the "Emergency open" button in the stop execution window.	User action
User: Exit	The user has chosen the "Exit" option of the pop-up menu.	User action
User: Fixed tube selected	The user had selected to fix the tube diameter	User action
User: Manual tube ID	The user had identified a tube manually.	User action
User: Must select tube diameter	The user wants identify tubes without selecting a fixed tube diameter.	User action
User: Open	The user has clicked the "Open door" button in the status window.	User action
User: Open door	The door was opened due an user request.	User action
User: Quality Control	The user has chosen the "Quality Control" option of the pop-up menu.	User action
User: Query Error Cancel	The user has clicked the "cancel" button of a error red window.	User action
User: Query Error Cancel all	The user has clicked the "cancel all" button of a error red window.	User action
User: Query Error Manual replacement	The user has clicked the "manual replacement" button of a error red window.	User action
User: Query Error Retry	The user has clicked the "retry" button of a error red window.	User action
User: Query Error Stop	The user has clicked the "stop" button of a error red window.	User action
User: RETRY TestWiz	The user has pushed the arrow button in the selection test window to retry the samples, reagents or cards identification.	User action
User: START	The user has clicked the arrow button in the status window.	User action
User: START , Cards Incidences ...	The user has clicked the question mark button in the status window when there was cards incidences.	User action
User: START , Reagents Incidences ...	The user has clicked the question mark button in the status window when there was reagents incidences.	User action
User: START , Sample Incidences ...	The user has clicked the question mark button in the status window when there was sample incidences.	User action
User: START TestWiz	The user has pushed the arrow button in the selection test window to begin a new execution.	User action
User: START Urgency	The user has pushed the arrow button in the selection test window to begin a urgency execution.	User action
User: Stop all	The user pushed the "stop!!" button in the stop execution window.	User action
User: Stop pipeting	The user pushed the "stop" button in the stop execution window.	User action
User: Stop pipeting and incubation	The user pushed the "stop!" button in the stop execution window.	User action
User: Stop request	The user request to stop	User action
User: Stop tubes identification	The user pushed the "stop" button in the selection test window to stop the samples identification started with the "Identify" button.	User action
User: Tubes Identification	The user has pushed the "Identify" button in the test selection window.	User action
Vial detected in position %d	A vial was detected in the shown position of the reagents carousel	Execution trail
Waiting for Door Closed	The door must be closed and are the execution is waiting for the user close it.	Execution trail
Waiting for timer in batch %d, subbatch %d	The execution stops waiting for a timer.	Execution trail
Washing block	A washing block begins.	Execution trail
Well %s	A well is dispensed.	Execution trail
Worksheet launched by user %s	User who begins the worksheet.	Execution trail
X cam%d : %ld	X calibration of a camera.	Execution trail
X Calibration for centrifuge: %d	Logged when a the centrifuge was calibrated	Execution trail
X Centri : %6ld	X calibration of the gripper in the centrifuge.	Execution trail
X Incubator: %6ld	X calibration of the gripper in the incubator.	Execution trail
XerErrCod %d	An error code was returned by the "Cavro" plunger. See the technical manual for a error description.	Errors trail
XYZErrCod %s	An error code was returned by the "XYZ" card. See the technical manual for a error description.	Errors trail
Y cam%d : %ld	Y calibration of a camera.	Execution trail
Y Calibration for centrifuge: %d	Logged when a the centrifuge was calibrated	Execution trail
Y Centri : %6ld	Y calibration of the gripper in the centrifuge.	Execution trail
Y Incubator: %6ld	Y calibration of the gripper in the incubator.	Execution trail
Z Calibration for centrifuge: %d	Logged when a the centrifuge was calibrated	Execution trail
Z Centri : %6ld	Z calibration of the gripper in the centrifuge.	Execution trail
Z Delay Centri : %d	Gripper delay for the centrifuge.	Execution trail
Z Delay Incub : %d	Gripper delay for the incubator.	Execution trail
Z Incubator: %6ld	Z calibration of the gripper in the incubator.	Execution trail
Z Needle : %6ld	Z calibration of the needle in the incubator	Execution trail
Z Reader : %6ld	Z calibration of the gripper in the reader.	Execution trail

GRIFOLS Diagnostic Grifols, S.A.	Reference:	REGD-0000847	Version:	3.0, CURRENT
	Status:	Effective	Effective Date:	26-May-2011
Title:	SM DG57-DG517. Log file description.			Page 9 of 9

DOCUMENT CHANGE HISTORY

Document change history			
Version	Date	Change	Rationale
2.0	2005/05/09	<i>Document REGD-0000847_v2 replaces document 57077V0M000.</i>	Document reference changed.
3.0	See document header	Added `previous version `data`, `change` and `rationale` inside document change history section.	Update change history for previous document history modification.
3.0	See document header	<i>Document REGD-0000847_v3 replaces document REGD-0000847_v2.</i>	Document version changed.
3.0	See document header	Removed software version from log example.	Document updated for current and following versions.
3.0	See document header	Document title.	Commercial reference substituted by project reference.
3.0	See document header	Header title.	Increase information to easy document search.