



Ppm Checklist of Chemistry Analyzer

Work Order		Ppm Date	
Analyzer Model		Hospital	
Serial Number		Department	
Windows OS		SW Version	

Replace parts

No.	Action Taken	Check	Comment
1	Replace Lamp		
2	Replace Reagent/Sample probe		
3	Replace Mixer Bar		
4	Replace Syringe Plunger		
5	Replace DI water filter		

Adjustment/Check

No.	Action Taken	Check	Comment
1	Sample/Reagent probe adjustment		
2	Mixer Bar adjustment		
3	Check cooling		
4	Check Heating		

Lubrication/Cleaning

No.	Action Taken	Check	Comment
1	Gas spring lubrication		
2	Dispenser rod lubrication		
3	Clean Reagent/sample carousel		
4	Clean surface of analyzer		
5	Clean wash station		

Testing

No.	Action Taken	Check	Comment
1	Initialize machine and prime		
2	Run samples to check machine		

Engineer Name

Signature