



Luminix PM CHECK LIST

No.	Procedure	6 M PM	12 M PM
1	Assy, Sheath Filter Tubing change	☐	☐
2	Syringe, 500µL, Cavo change	☐	☐
3	Tubing, PEEK, 1/16 OD, .010 ID, 10" long change	☐	☐
4	Assy, Tubing, Common Valve to Pump, 50.0" PTFE change	☐	☐
5	Assy, Tubing, Outside Sample Tube, 5.0" PTFE	☐	☐
6	Filter, Air, HEPA, 1-5 micron, 3/16" barb out change	☐	☐
7	Filter, Fan, XY, Foam Media, 30 PP change	☐	☐
8	Operating voltage check	☐	☐
9	Heating check	☐	☐
10	Cooling check	☐	☐
11	Washing station cleaning	☐	☐
12	Mechanics cleaning / lubrication	☐	☐

Date:

Place:

Service done by:
