

AmiCORE Installation/Preventive Maintenance Checklist

AmiCORE S/N: _____	Date: _____
Customer: _____	SW Version: _____
Service #: _____	FSE Name: _____

Instructions:

1. For detailed information, including required test points, refer to the AmiCORE Service Manual.
2. For installations, complete Sections I and III. For Preventative Maintenance, complete Sections II and III.

Check box note: P = Pass, F = Fail, N/A = Not Applicable

I. Installation Only	NA ☐	
3.2 Instrument Unpacking & Assembly	P ☐ F ☐	3.14 Interface Verification
6.1 Touch Screen Calibration	P ☐ F ☐	3.15 Weight Scale Verification
6.2 Weight Scale Verification & Calibration	P ☐ F ☐	3.16 Leak Detector Verification
6.3 Centrifuge Pressure Sensor Calibration	P ☐ F ☐	3.17 Temperature Sensor Verification
6.4 Interface Detector Assembly Alignment Calibration	P ☐ F ☐	3.18 Cuff Verification
		3.19 USB Verification
		3.20 Power Verification
II. Preventative Maintenance Only	NA ☒	
4.2 Parts Replacement P.M.	P ☐ F ☐	3.21 Battery Verification
		3.22 Simulated AmiCORE Procedure
III. Installation & Preventive Maintenance		
3.3 Visual Inspection	Completed ☒	
3.4 Power UP System Test	P ☒ F ☐	
3.5 Enable Service Mode	P ☒ F ☐	
3.6 Pump Verification	P ☒ F ☐	
3.7 Cassette Motor Verification	P ☒ F ☐	
3.8 Clamp Verification	P ☒ F ☐	
3.9 Valve Verification	P ☒ F ☐	
3.10 Pressure Sensor Verification	P ☒ F ☐	
3.11 Air Detector Verification	P ☒ F ☐	
3.12 Door Lock Verification	P ☒ F ☐	
3.13 Centrifuge Verification	P ☒ F ☐	

Comments: _____

Signature: _____ Date: _____