



**Compo Seal**  
Yearly MAINTENANCE  
CHECKLIST



|               |  |            |  |
|---------------|--|------------|--|
| Work Order No |  | PPM Date   |  |
| Region        |  | Hospital   |  |
| Serial No     |  | Department |  |

| Action                                |                                                                                                                                  | Checked / Done           | Remarks |
|---------------------------------------|----------------------------------------------------------------------------------------------------------------------------------|--------------------------|---------|
| <b>General checks</b>                 |                                                                                                                                  |                          |         |
| 1                                     | Visual inspection: instrument general status.                                                                                    | <input type="checkbox"/> |         |
| 2                                     | Mains Inlet Status                                                                                                               | <input type="checkbox"/> |         |
| <b>Decontamination &amp; Cleaning</b> |                                                                                                                                  |                          |         |
| 1                                     | Turn off the instrument, disconnect the power cord.                                                                              | <input type="checkbox"/> |         |
| 2                                     | If there is spillage, wipe it with disposable absorbent material (for example, paper towels, gauze, or cloth).                   | <input type="checkbox"/> |         |
| 3                                     | Wipe the entire instrument with disposable towels or cotton balls dipped in a 0.5% (v/v) aqueous solution of sodium hypochlorite | <input type="checkbox"/> |         |
| 4                                     | Rinse the surface with a soft detergent and water to eliminate odors and any residual solution.                                  | <input type="checkbox"/> |         |
| 5                                     | Dry the surface.                                                                                                                 | <input type="checkbox"/> |         |
| 6                                     | Dispose of contaminated materials used during decontamination in a container for biological waste.                               | <input type="checkbox"/> |         |

|                  |  |
|------------------|--|
| Service Engineer |  |
|------------------|--|