

URISED UNINSTALLATION CHECK LIST

Customer Name : _____ Inst SN# _____

S.NO.	Description	Status	Comment
1	Remove the unused cuvettes		
2	Empty and clean The used Cuvettes Bin		
3	Prime the instrument with air		
4	Empty and clean the waste and Water bottles		
5	Remove all the accessories and computer		
6	Clean and disinfect the instrument		
7	Lock the robot arm		
8	Pack the accessories		
9	Fill the packing list		

Checked By: _____ Signature: _____ Date: _____