

GRIFOLS Diagnostic Grifols, S.A.	Reference:	REGD-0000774		
	Version:	5.0; Most-Recent; Effective; CURRENT		
	Status:	Effective	Purpose:	Uncontrolled Copy
Class:	Control del Diseño		Effective date:	28-Apr-2016
Title:	SM DG57. Technical specifications			

DG57 Service Manual

TECHNICAL SPECIFICATIONS

Document description
This document describes the technical specifications of the instrument. That includes: - Operating and environmental. - Computer requirements. - UPS recommendation. <i>Note: Document REGD-0000774_v5 replaces document REGD-0000774_v4.</i>

GRIFOLS Diagnostic Grifols, S.A.	Reference:	REGD-0000774		
	Version:	5.0; Most-Recent; Effective; CURRENT		
	Status:	Effective	Purpose:	Uncontrolled Copy
Class:	Control del Diseño		Effective date:	28-Apr-2016
Title:	SM DG57. Technical specifications			

Specifications				
Operation and Environmental specifications				
1.				
	MODEL:	DG-57 (PC not included)		
	SUPPLY:	Voltage:	100-120/220-240 V~	
		Frequency:	50-60 Hz	
		Input power:	300 W	
		Fuses:	2 x T4A L 250V, 5x20 mm	
	PROTECTION AGAINST ELECTRIC SHOCK:	Class:	I	
	INSTALLATION CATEGORY:	Overvoltage category II (local level, appliances, portable equipment, etc.)		
	DIMENSIONS (mm):	600 (depth) x 1000 (width) x 650 (height, with door closed)		
	WEIGHT (Kg):	Approximately 87 kg		
	MATERIALS:	Packaging: wood, ethafoam Equipment: aluminium, stainless steel, iron, polyurethane (expanded foam), glass, electronic components (silicon, copper...).		
	OPERATING CONDITIONS:	Indoor use		
		Temperature:	15°C to 28°C	
		Maximum relative humidity non condensing:	80%.	
		Maximum altitude:	2000m.	
		Mains supply maximum voltage fluctuations:	±10% of nominal voltage.	
	TRANSPORT AND STORAGE CONDITIONS:	Temperature:	-10°C to 60°C	
		Maximum relative humidity non condensing:	85%	
	CENTRIFUGE:	Maximum load:	12 cards	
	WORKING PRESSURE OF BOTTLES:	approx. 180 mbar (pressure) < -300 mbar (vacuum)		

GRIFOLS Diagnostic Grifols, S.A.	Reference:	REGD-0000774		
	Version:	5.0; Most-Recent; Effective; CURRENT		
	Status:	Effective	Purpose:	Uncontrolled Copy
Class:	Control del Diseño		Effective date:	28-Apr-2016
Title:	SM DG57. Technical specifications			

Specifications

Computer requirements

2. In case of using the Matrox 4Sight GP computer, the default configuration is:

- a. Operating system: Windows Embedded Standard 7.
- b. RAM memory: 4GB DDR3
- c. Hard disc: 500GB SATA
- d. 2 USB 2.0 ports (front)
- e. 4 USB 3.0 ports (rear)
- f. 3 Serial ports (including 2 port PCI Express board)
- g. Matrox Morphis eDual standard (PCI Express) included

Note: A DVI-I to VGA adaptor is included in the PC package in order to be able to connect it to VGA screens.

3. Those are the minimum requirements for a standard PC to work with the software.

- a. Processor: Pentium II or higher
- b. Memory: minimum 64Mb
- c. Operating system: Windows 2000, XP and Windows 7 Professional 32 bits
- d. Other:
 - i. 1 free PCI or PCI Express slot, for Matrox Meteor II or Matrox Morphis image acquisition board (respectively).
 - ii. Printer port.
 - iii. At least 1 COM port.

UPS recommendations

4. In case of using a UPS unit, ensure it complies these specifications:

- The instrument is designed to support mains microcuts. Thus, it is not necessary to install a UPS unit only for this purpose.
- Output voltage: compatible with the one of the equipment.
- Power: a 800VA UPS is enough to give power to the instrument and the computer with monitor attached to it. The power consumption of the instrument is 150VA.
- Use always a certified UPS (Electromagnetic compatibility and Electric safety regulations).

GRIFOLS Diagnostic Grifols, S.A.	Reference:	REGD-0000774		
	Version:	5.0; Most-Recent; Effective; CURRENT		
	Status:	Effective	Purpose:	Uncontrolled Copy
Class:	Control del Diseño		Effective date:	28-Apr-2016
Title:	SM DG57. Technical specifications			

DOCUMENT CHANGE HISTORY

Document change history			
Version	Date	Change	Rationale
2.0	2005/05/09	Updated operational and environmental specifications.	Operational and environmental specifications outdated
2.0	2005/05/09	Document REGD-0000774_v2 replaces document 57002V0M000.	Document's reference changed.
3.0	May 2011	Replaced commercial name by project name.	Delete commercial name.
3.0	May 2011	Document REGD-0000774_v3 replaces document REGD-0000774_v2.	Document version changed.
3.0	May 2011	Minimum relative humidity non condensing updated.	Operational and environmental specifications updated.
3.0	May 2011	Computer requirements updated.	Operating system updated.
3.0	May 2011	Header title has changed.	Previous header title has been substituted by project and document title.
4.0	December 2014	Computer requirements	New PC 4Sight GP
5.0	April 2016	Matrox 4Sight GP default configuration updated..	Technical specifications update
5.0	April 2016	Operating conditions updated	Technical specifications update